

Product datasheet for **MR213986L3V**

Strip2 (NM_001037740) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Strip2 (NM_001037740) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Strip2
Synonyms:	D330017J20Rik; Fam40b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001037740
ORF Size:	2448 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR213986).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001037740.1 , NP_001032829.1
RefSeq Size:	5286 bp
RefSeq ORF:	2451 bp
Locus ID:	320609
UniProt ID:	Q8C9H6
Cytogenetics:	6 A3.3
Gene Summary:	Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)