

Product datasheet for **MR214014**

Pih1h3b (NM_177921) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pih1h3b (NM_177921) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pih1h3b
Synonyms: E230019M04Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR214014 representing NM_177921
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAATCGGAAATTGCAGAAGCTGAGAAGGTCAAACAGGGAATTTGGAAGCAAAAAGCATGGAAGCTA
AAAATATGCAATTTGAGATGCTGTCTTCTGCCTCAACCTGCAAGCACTGGCCAGTCTCCTGAATCCTGA
AGAAGAGGATGATTTTGATTTTGAACAGGGACAATGTACATCTACTATTAGAGCTATGGGCCCTGGGAAT
ATTGGACCGCCTAAAGCAAAAGAGTCAAAGCCATTCTGAACCCAGGAGTGATGAAAGTAAAACATCT
GGAATCCAGAAGAGGTTCCAGAAGGAGGAGAGCATGATGATATATGGGATGTTTCGAGAAATCCCAGATTA
TGAGATTGATTCCAACAACTGTGGAACTGAGGACATATACCTAGGATTGACAGGAAAAGGATCCTTCA
ACAGCATGCTGCCAGGAGCTGGTGGTGAATAAATACTGCCAATACGAACCTTCTGAAATTCAAATTG
ACATCCAGGAAATGCTCCTTGATCTTCGCACTCCTAAAAAGAAGCTGTTAGTAACTTTCCCCAGCCTGT
GGAACGGAACAGTGCTAGAGCATCCTATATCTGGGAAGCTGAAACACTTGAAGTCAGGATGACCCCTCAA
AGGGATTTGGATTTAATATTTC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR214014 representing NM_177921
Red=Cloning site Green=Tags(s)

MESEIAEAEKVKTGNLEAKSMEAKNMQFEMLSSASTLQALASLLNPEEEDDFDFEQGQCTSTIRAMGPGN
 IGPPKAKESKAIPEPRSDESENIWNPEEVPEGGEHDDIWDVREIPDYEIVFQQTVGTEDIYLGITGKDPS
 TACCQELVVKIKLPNTNPSEIQIDIQEMLLDLRTPKKLLVNFPPVVERNSARASYIWEAETLEVRMTLQ
 RDLDFNIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_177921

ORF Size: 654 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_177921.3](#), [NP_808589.1](#)

RefSeq Size: 1509 bp


RefSeq ORF: 657 bp

Locus ID: 331537

Cytogenetics: X F1

MW: 24.5 kDa

Product images:



Circular map for MR214014