

Product datasheet for **MR214066L4V**

Ino80d (NM_001081436) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ino80d (NM_001081436) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ino80d
Synonyms:	7330405I11; A430093A21Rik; BC010584; OTTMUSG00000002594
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001081436
ORF Size:	2748 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214066).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081436.2 , NP_001074905.1
RefSeq Size:	13438 bp
RefSeq ORF:	2751 bp
Locus ID:	227195
UniProt ID:	Q66JY2
Cytogenetics:	1 C2
Gene Summary:	Putative regulatory component of the chromatin remodeling INO80 complex which is involved in transcriptional regulation, DNA replication and probably DNA repair. [UniProtKB/Swiss-Prot Function]


[View online »](#)