

Product datasheet for **MR214079**

Gtpbp8 (NM_025332) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Gtpbp8
Synonyms:	0610037H22Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR214079 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCAGCGGCTAGGTTGAGTCATAGAATGGGAAGGCTCCTTGAGAAAGCGCCAGCACTGGGGCCTTGGACCCGGGTTTACAGCACGTCCCCTGCCTTCGAGAAGTGCTAAGGCTGCCACAGAAACAGCTAACGAAGGTAGTGTACCCGCTGCGGGAACGGGGCAGCACCTCGCCGCGGACTCGGGGCTGGCCTTATCGAACAGAGACTCTTCGACCCGAGCCTGGAGGACATCGGGCGGGCGGAGAGCATCTTTGAGGCCACCGCCCGGAACCGCATCGAGTACCTCAGCTCCGCCGTCCGCCTCGACCACGCCCCGAGTTTGACGACGCGGAGGTGTGTTTCATAGGCCGAAGCAATGTTGAAAAATCCTCTTGATAAAGGCCTTATTTCCCTGGCCCCTGATGTAGAAGTCAGGATCTCAAAAAACCGGGGCACACAAAGAAAATGAATTTTTTCAAAGTTGGAAAACATTTTACATTGGTGGACATGCCAGGGTATGGCTATAGAGCACCCGAAGACTTCGTTGACATGGTAGAGACTTACCTAAAGGACGAAATAATTTGAAGAGGACATTTTATTAGTGGATAGTGTGTTGGAATTACAAAATTAGACAATATTGCCATAGAAATGTGTGAAGAATTTGCATTGCCTTATGTGATGATATTACAAAAATCGACAAATCTTCCAGGGGTATCTTTTAAACAAGTGCTCCAGATTCAGAAGTTTGTTAACACGCAGACACAAGGCTGTTTTCACAGCTGTTTCTATAAGTGCAGTGACCAATTCGGAGTCCATCTGTTGAAATGCTTTATCGCAGATATAACAGGAAGTCTTAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR214079 protein sequence
Red=Cloning site Green=Tags(s)

MAAARLSHRMGRLLKAPALGPWTRVYSTSPFAEVLRLPQKQLTKVVYPLRELQGHAAADSGPLIEQR
LFDPSLEDIGRAESIFEATARNRIEYLSAVRLDHAPSLQQPEVCFIGRSNVGKSSLIKALFSLAPDVEV
RISKPGHTKKMNFVKGHFTLVDMPGYGYRAPEDFVDMVETYLKERNNLKRTFLLVDSVVGITKLDNI
AIEMSEEFALPYVMILTKIDKSSKGYLLKQVLQIQKFVNTQTQGCFFPQLFPIASVNTSGVHLLKCFIADI
TGSLCK



TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN:

NM 025332

ORF Size:

855 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_025332.3](#), [NP_079608.1](#)

RefSeq Size: 1318 bp

RefSeq ORF: 858 bp


Locus ID: 66067

UniProt ID: [Q9CY28](#)

Cytogenetics: 16 B4

MW: 31.9 kDa

Product images:



Circular map for MR214079