

Product datasheet for **MR214080**

Ugt2b36 (NM_001029867) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Ugt2b36
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >MR214080 representing NM_001029867
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCTTTGGAAGTGGATTCTGCTCTGCTGCTGCTGCAGATAAGTTGCTGCTTCAGATCTGCAAAGTGTG
 GGAAGGTGTTGGTATGGCCGGTGGACTACAGTCACTGGATGAATATAAGATAATACTGGATGAACCTAA
 ACAGAAGGGCCATGAAGTCACTGTTTTGAGACCTTCAACTTCCATTTTTCTTGATCCCAAAAAATCGCCT
 GGCTTAAGTTGAACTTTTCTACATCTTTCAGTAATGATGTCATGGAATATTTTTGCAAAGGCTG
 TGGAAAGATGGACTTATGAGGTTCCAAGAGATACATGTTTATCATATTCTCTTTATTACAAAATATATT
 TGATGAATATTCTGATTACTGTTTAACTCTTTGCAAAGACACAGTTTCAAACAAACAGCTCATGGCAAAA
 CTCCAGGAATCCAAGTTGATGTCATTTTGTGAGATGCTATTGGTCCCTGTGGGAGCTCATAGCTGAAC
 TGCTCCAAATTCTTTTCTATATAGTCTTCTGCTCTCTCTGGCTACTATCTTGAAAAATACAGTGGAGG
 ACTCCCTCTCCCTCCTTATGTACCTGTAATTTGTCAGGATTAAGTGGTCAAATGACATTCAAGAG
 AGGGTAAAGAATATGATATGTATGCTTTATTTGACTTTTGGTTCCAGACATTTAGGGAGAAGAAATGGG
 ATCAGTTTTACAGCGAACTTTGGGAAGGCCAACCTTAATTGAGACAATGGGCAAGCAGAAATGTG
 GCTCATTAGATCCTATTGGGATTGGAGTTTCTCACCCAACCTTACCAATGTTTACTATGTTGGAGGA
 GTCCATTGCAAACTGCTAAACCTTGCTAAGGAAATGGAAGACTTTGTTGAGAGCTCTGGAGAGCATG
 GTGTGGTGGTGTCTCTGCGGTCAATGGTCAGTAATATGACAGAAGAAAAGGCCAATGCAATTGCATG
 GGCTCTTGCCAGATTCACAAAAGGTTCTTTGAGATTGACGGCAAACTCCAGCCACCTTAGGACCC
 AATACCAGAATCTACAAGTGGCTTCCCAAAATGATCTTCTGGTCATCCAAAACCAAGCCTTTATAA
 CTCATGGTGGCGGAATGGACTCTACGAGGCGATCCATCATGGAATCCCTATGATTGGCATTCTTTGTT
 TGGAGAACACATGATAACATTGCTCACATGGTGGCCAAAGGAGCAGCTGTACATTGAATATCAGAACG
 ATGTCAAGGTCAGATTTGCTCAATGCACTGGAGGAAGTCATAGACAATCCTTTCTATAAAGAGAATGCTA
 TGTGGTTGTCAACCATTCACCATGACCAGCCTATGAAGCCCTGGACAGAGCAGTCTTCTGGATTGAGTT
 TGTCTGCGCCACAAAGGGGCAAGCATCTGAGACCACTTGATACAACCTCACCTGGTACCAGTACCAC
 TCTCTGGATGTGATTGGATTCTGCTTGTGTTGACATTATTGTAGCCCTCATTGTAAAGTGCTTCT
 TGTTTGTACCATTCTTTGTAAAGAAAGAAAAGAAATGAAGAATGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR214080 representing NM_001029867
 Red=Cloning site Green=Tags(s)


MLWKWISALLLQISCCFRSAKCGKLVLPVDYSHWMNIKIILDELKQKGHEVTVLRPSTSIFLDPKKSP
 GLKFETFTSFSNDVMEIIFAKAVERWTYEVPRDTCLSYSPLLQNIFFDEYSYCLTLCKDVTNSKQLMAK
 LQESKFDVILSDAIGPCGELIAELLQIPFLYSLRFSFGYYLEKYSGLPLPPSYVPVILSGLSGQMTFKE
 RVKNMICMLYDFWFQTFREKKWDQFYSETLGRPTTLIETMGKAEMWLIIRSYWDLEFPHPPTLPNVYYVGG
 VHCKPAKPLPKEMEDFVQSSGEHGVVVFSLGSMVSNMTEEKANAIWALAQIPQKVLWRFDGKTPATLGP
 NTRIYKWLQNDDLGHPTKAFITHGGANGLYEAIIHHGIPMIGIPLFGEQHDNIAHMAVAKGAAVTLNIRT
 MSRSDLLNALEEVIDNPFYKENAMWLSTIHHQPMKPLDRAVFWIEFVMRHKGAKHLRPLAYNLTWYQYH
 SLDVIGFLLAFVTFIVALIVKCLFVYRFFVKKEKMKNE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9015_g07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001029867

ORF Size: 1590 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001029867.1](#), [NP_001025038.1](#)

RefSeq Size: 1888 bp

RefSeq ORF: 1593 bp
Locus ID: 231396
Cytogenetics: 5 E1
MW: 61.5 kDa

Product images:



Circular map for MR214080