

## Product datasheet for **MR214084**

### Sult3a2 (NM\_001101452) Mouse Tagged ORF Clone

#### Product data:

|                           |                        |
|---------------------------|------------------------|
| Product Type:             | Expression Plasmids    |
| Tag:                      | Myc-DDK                |
| Symbol:                   | Sult3a2                |
| Synonyms:                 | EG215895; Gm4794       |
| Mammalian Cell Selection: | Neomycin               |
| Vector:                   | pCMV6-Entry (PS100001) |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |

**ORF Nucleotide Sequence:** >MR214084 representing NM\_001101452  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGACAATAAAGATGAAGATATTTTTAAATTCAAAGGCTATAATTTCTCAGAGCGTTAGTTAGTATGG  
ATATAGTGAAAAATATGGAAGATTATGAAATTCGGGATGATGACATCTTATAGTCACATACCCAAAGTC  
TGGTACCACCTGGACTCAACAGATCCTCTGCTTAATTTGTTTTGAGAGTCACCGGAATGGAAGTAAAAAC  
ATCGCAACAATAGATAGAATACCATCCTTTGAGTACAATATTCACAAATTAGACTATGCCAAATGCCAT  
CCCCTCGAATCTTCACTTCCACATTCCATATTACTTAGTACCAAAAGGTCTCAAGGACAAAAAGCCAA  
AATTATTTATGTCTACAGAAATCCTAAAGATGTTTTGATCTCCTTTTTTCATTTCTCAAATTTGATGGTT  
GTACTCAAAGCTTCAGATACTCTGAAAAATTTATGCAAAGATTTCTAAATGGAACTTGGTAGGAAGCC  
TTTGGTTTGACCACATAAGAGGCTGGTATGAACACAGACATGATTTCAATATAATGTTTCATGAGCTTTGA  
AGACATGAAGAAGGACCTCAGAAGCTCAGTGCTTAAATCTGCAGTTTTCTGGAGAAAGAACTGAGTGAA  
GAAGATATGGATGCTGTTGTGAGACAAGCTACCTTTGAGAAAAAGAACTGACCCACGAGCAAAACATG  
AACACATTATAAAGATGAAGTACGGAACAAGAAATAACGGGAGCTTTCTTCGCAAGGTATCGTTGG  
AGACTGGAAACATTATTTGACTGTGGATCAAAGTGAAGATTTGACAAGATTTTCATATGAACATGAAA  
AACATCCCCTTAAATTCATTTGGAATATAAATGAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR214084 representing NM\_001101452

Red=Cloning site Green=Tags(s)

MDNKDEDIFKFKGYNFLRALVSM DIVKNMEDYEIRDDIFIVTYPKSGTTWTQQILCLICFESHNGTEN  
IATIDRIPSF EYNIHKLDYAKMPSPRIFTSHIPYYLVPKGLKDKKAKIIYVYRNPKDVLISFFHFSNLMV  
VLKASDTLENFMQRFLNGNLVGSLWFDHIRGWYEHHRDFNIMFMSFEDMKDLRSSLVKICSFLEKELSE  
EDMDAVVRQATFEKMKADPRANNEHIKDELGTRNNTGAFLRKGI VGDWKHYLTVDQSERFDKIFHMMNK  
NIPLKFIWNINEK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

[https://cdn.origene.com/chromatograms/mm9081\\_d01.zip](https://cdn.origene.com/chromatograms/mm9081_d01.zip)

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM\_001101452

ORF Size:

879 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_001101452.1](#), [NP\\_001094922.1](#)

**RefSeq Size:** 882 bp


**RefSeq ORF:** 882 bp

**Locus ID:** 215895

**Cytogenetics:** 10 B1

**MW:** 34.8 kDa

**Product images:**



Circular map for MR214084