

Product datasheet for **MR214085**

Olfr1153 (NM_146640) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Olfr1153
Synonyms:	MOR177-7
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR214085 representing NM_146640
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCAGTTTATAAAACAAAGAAGAATGGAGGTTGGAACTGTTCAGCCACTGAGTTTCTTCTTTGGGAA
TTACCAATAATCCTGTGATCAAGGTGATTTTATTACTACTTTTCTAATTGTTTATCTCATAATCTTAT
TGAAAATCTGGGATGATCATTTTGATAAGAATGAACCTCCAGCTTCACACACCGATGTACTTTTTCTC
AGCCACCTCTCTTCTCTGATATCTGTTATTCTACTGCAGTTGGACCAAGATGTTGGTGGCCTTATAT
TCAAGAACAAATTCATTCCATTATTGACTGTGCTGTGCAGTTCTTTATCTTCTGCATTTTACCGATGC
TGAGTGTGTGCTACTGGCAGTGATGGCTTTGACAGGTATAAGGCCATCAGCAACCCCTGATGTATGCG
GTAGATATGTCCAGAGTGGTATGCTATCAGCTCCTGGCTGTGTTTATCTGGTTGGAATGGTAGATGCTT
TGACACATAACAAGTTGACATTTCACTTGTGTTTCTGTCAGTCAAAGAAATCAACCATTCTTCTGTGA
TGTTCTCTCTCTCTATTACTATCTTGTTTCAGATACAGAGGTCAATGAAGTATGTTATTTACCCCTATTT
GGGTTTATTGAGCTGAGCACCATCTCGGGAGTCTGGTCTCTTACTGTTACATCATCTCATGTTTGA
AGATTGCTCTGCTGAGGGCAGATTCAAAGCATCTCCACCTGTACCTCCCATCTAAGTGTGTTGCAAT
TTTCCAGGGAACAATGCTCTTCATGATTTCCGGCCAAGTTCTGCTTATCCCTGGATCAAGACAAATG
ACCTCATTATTCTATACTCTTGATCCCATGCTCAACCCGTTGATTTACAGCCTCAGGAATAAGGATG
TAAAGCAAGTGTAAGGTCCTGAAAAGTAGAATTTTATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

```
>MR214085 representing NM_146640
```

Red=Cloning site Green=Tags(s)

MQFIKQRRMEVGNCSEATEFLLLGITNPNVIVILFTFLIVYLIIL IENLGMIIIRMNSQLHTPMYFFL
SHLSFSDICYSTAVGPKMLVGLIFKNNSIPFIDCAVQFFIFCIFTDAECVLLAVMAFDTRYKAISNPLMYA
VDMRSRVVCYQLLAVVYLVGMVDALHTTTLTHFLCFCQSKENINHFFCDVPPLLLLSCSDTEVNELVIFTLF
GFIELTSISGVLVSYCIISVLLKIRSAEGRFAKFCSTSHLTVAVIFQGTMLMFPYFRPSSAYSLLDQDKM
TSLFTYIVIMPLNPLYISLRNDKVCASVKRSLKSRISF

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Chromatograms:



https://cdn.origene.com/chromatograms/mm9084_c12.zip

Restriction Sites:

Sqfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN:

NM_146640

ORF Size:

954 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_146640.2, NP_666851.2</u>
RefSeq Size:	957 bp
RefSeq ORF:	957 bp
Locus ID:	258633
Cytogenetics:	2 D
MW:	36.1 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR214085