

Product datasheet for MR214142

Prr29 (NM_029345) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prr29 (NM_029345) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Prr29
Synonyms: 2310007L24Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR214142 representing NM_029345
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCTCCCTCGGGAGTGGAACTGGAGTGGTGTCCCAACACAGAGCACAGCCCTATGCCCTGGGTGA
 CCATCTTGACGCCCTTCCTTGGACAGTCCCATCTTCTCAGCCTCAGCATAACCGGGTCAAGGAAGACCT
 GCTGGAGCTGATGATGCTACAGAACGCACAGATGCACCAGCTGCTACTGGCCAACTGGTGGCAGATGCC
 CTGAACCCAGGGCCAGAGTGGCCAGCCACCGGTCCACACCGACAGTCACGAGGAGCAGATGGAGGAAG
 AGATGGAGATGCAGGAGCAAGAGCCTTTGGTTTACCACCATCACTACTTGCCTTGCCAGTGACTCCCTT
 GGGTCCCATGTCCTCTGGCCAGCCTCCTTCTCCCTGTTCCCCACATCAGCCTCCCTGGCAGGGTGAA
 CCGAGAATTCAGCACCAGCCCCCTGCCTCCAGGCAAGAGGAGGTGAGAGATGTGCCCCACCCCCACCCC
 CCAGTGAACGGGGACTGTTGGTGCAGATGTACCTCCAGCTTCAGACTACTATGATGCGGAGAGCCTGCC
 A

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR214142 representing NM_029345
 Red=Cloning site Green=Tags(s)

MASLGSNWSGVPTQSTAPMPWVTILQFPFWTPSSQPQHNRVKEDLLELMMLQNAQMHQLLLGQLVADA
 LNPGPEWSPPVHTDSHEEQMEEEMEQEPLVYHHHYLPCPVTPPLGPMASLWPASFLPVPPHPQPPWQGE
 PRIQHQPASRQEVRDVPVPPPPSATGTVGADVPPASDYDAESLP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/mm9079_f05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_029345

ORF Size: 561 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029345.1](#), [NP_083621.1](#)


RefSeq Size: 789 bp

RefSeq ORF: 564 bp

Locus ID: 75573

UniProt ID: [B1ARI9](#)
 Cytogenetics: 11 E1
 MW: 20.8 kDa

Product images:



Circular map for MR214142