

Product datasheet for **MR214271**

Olfr948 (NM_001011756) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr948 (NM_001011756) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr948
Synonyms:	MOR171-45; MOR171-51
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR214271 representing NM_001011756 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCAGCTGCAAACCATTCATAGTACTGAGTTCTTCTGGCTGGACTCTCAGAGAATCAAAAAGTTC
 AGCTGCCCTTTTCTCCTATTTGTAGCAATCTATCTGATCACTGTGGCAGGAAACCTGGGGATGATTGC
 ATTGATTGGGATCAGTTCCCACTTGCACACCCCATGTACTATTCCTCAGCAGTCTGTCCTTCATTGAC
 TTCTGTCAGTCTACAGTTGTTACCCCTAAAATGTTGGTGAGCCTTCTGACAAAGAAGAATCATCTCCT
 ACTCTGGATGCATGGTTCAGCTCTACTTCTTCATCAGTTTGGAACTGCAGAGTGCTACACATTAGCTGT
 AATGGCATATGACCGCTATGTTGCCATTTGTAACCCCTACGTTACAATGTAACCATGTCCTATCAGATT
 TACAGTTCCTTGATTTACAGGGGTATATTTATGCTGTGTTCTGTGCATCAGTAAACACCGGCTTTATAA
 TTAGGATTCAGTTCTGCAAGTTAAATGTGATCAACCACTATTTCTGTGATCTTCTCCCTGTTGAACT
 TGCATGCTCTAATACCTATATCAATGAAATATTGATTCTATTTTTGGTTCAGTAAACATCTGTGTCCCA
 ATGCTGACTGTTATTACTTCTACATCTTCATCATTGCCAGCATCCTCCGATTCGCTCCAGTGAGGGCA
 AGTTCAAAGCCTTCAGTACTGCAGTCCACATCTCTGCTGTTGCTATCTTGATGTTCTACTGCATT
 CACATACTTACAGCCTTCATCAGTGAGTTTGGTGGACCAAGGAAAGTGCCTCTGTGTTTATACTACT
 GTTGTACCTATGCTGAACCCATTGATCTACAGCTTGAGGAATAAGGATGTCACCTTGCAATTGAAAAGAA
 TTCTTGAACAAAAAGGCTTTATGTAAGCAGAAGTCGTGTGAAGATTATTAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR214271 representing NM_001011756
 Red=Cloning site Green=Tags(s)


MAAANHCIIVTEFFLAGLSENQKVQLPLFLLFVAIYLITVAGNLGMIALIGISSHLHTPMYYFLSSLSFID
 FCQSTVVTPKMLVSLTKKNIISYSGCMVQLYFFISFGTAECYTLAVMAYDRYVAICNPLRYNVTMSYQI
 YSSLISGVYIYAVFCASVNTGFIIRIQFCKLVINHYFCDLLPLLKLACSNTYINEILILFFGSVNICVP
 MLTVITSYIFIIASILIRSRSEKFKAFSTCSSHISAVAILYGSTAFTYLPSSVSLVDQKQVSSVFYTT
 VVPMNLPLIYSLRNKDVTLALKRILEQKRLYVSRSRVKIIN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9065_h03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001011756

ORF Size: 963 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001011756.1](#), [NP_001011756.1](#)

RefSeq Size: 966 bp

RefSeq ORF: 966 bp


Locus ID: 257912

Cytogenetics: 9 A5.1

MW: 36 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR214271