

Product datasheet for MR214300

1700067K01Rik (NM_183097) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	1700067K01Rik (NM_183097) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	1700067K01Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR214300 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTT**CAGCCCC**CACTGACCAGCTCCGTTACCTGCTCAGGAAGGCGGATGACTTTCAAAGCTACTTGC
TCTACAGCAGGGACCGAGTACAGAAGGAACAGCTAGCCAAGGCCATGCCACCTTCTTAACAATGTGTGA
ACCATACTTCTTATACCTTGAGGCTGCAGCAAGGAGTGTCCCTCCCATCTATGGAGCC**TCGAGGAGTTG**
GTCCGGAAGGGTTGCTGGAGATCTCCAGCAGCTAACCTGCGCTTGGAA**CAGTTAGTTCTCATGTATG**
CCTCCTTTGGGTT**CGTGGAC**CTGGAGGAAACCAACCCTTTGAGCATCTCCTGCTTTTT**TCGCGGGAGGTT**
CTCCATCAGCCCC**CCCAGATGTGTCCATCTTCCGATACTGCACACCAGCTGCCTATACTGCCAGTCAC**
TTCCCAAGATACCTCTATAAGAAGATGCCTGGAACCTGGAAACCACACAGAGGCCAGCAGCAAGGAG
CAGACTCCCATGTGGATTACTACTTCTTGTGTTATCGAGATACATGGGAGGATGCAGGCCAGGGACCAGC
CAATTCTGCCCCAGATCCAGAACTCTGGTCCATTGGCCGATGGATGCCCTAGGACCAGCAGAGGAT
GATCTTGACTCGTGGATTTTGTGCCACAGCCACCCGGGACTACCAGCAGCTACTAACCATCGGCTTCG
AGGAGCCGTCGCATGTGCTGGCCACCGACCTGCTGGTGCAGATCCTCCTGGCCAGGCAGCCCCCGCCG
GCCCC**CAGCGCTGCCGGCCCGCGCTTGGGCAGCCCCGGCATCC**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >MR214300 protein sequence
Red=Cloning site Green=Tags(s)

MFSPITDQLRYLLRKADDFQSYLLYSRDRVQKEQLAKAMPTFLTMCEPYFLYLEAAARSVPPIYGALQEL
 VRKGLLEISQQLTLRLEQLVLMYASFGFVDLEETNPLSISCFFCGRFSISPSHDVSI FRYCTPAAYTASH
 FPRYLKMKRWNL ETTTEASSQGADSHVDYF L CYRDTWEDAGQGPANSCPQIQKLWSIGRWMP L GPAED
 DLDSWILCPQPPGDYQQLLTIGFEESHVLATDLLLVQILLGQAGPARPPSAAGPAAWAAPAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_183097

ORF Size: 819 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_183097.2](#), [NP_898920.2](#)

RefSeq Size: 1084 bp

RefSeq ORF: 903 bp


Locus ID: 73453

UniProt ID: [G3UZ3](#)

Cytogenetics: 8 C2

MW: 30.7 kDa

Product images:



Circular map for MR214300