

Product datasheet for **MR214302A1V**

Mouse Ppp1r36 (NM_001163103) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	Ppp1r36
Synonyms:	Gm70
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR214302 representing NM_001163103 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGTCAAGAGTGAGGCCATGTTACCGTCCCAGAGTTTTACTCCAGGAGAAAGCAGTTTGTGGGCAAT
CCTCTACCGGCTGGACCAATGTGGATTACGACTAGGAATGTGGTATTGAAAGACGAAACAAAACTCT
GGAAATTTAAAGTTTCACACCTGCAATAGAACTCAAGGAGAAAGGCAAGAAAGGGAAGGCGGTTCACTTT
GCTGAGATGGACAGCGGAGCTTCGGAAGGTTGACAGATAAAAGGTTTGCTCAAGGGATGAGAAGGCGG
CTAAGGCTTTGGAGAAACGAAGTCAGCAGGGCAATGTCACGCTGGATGACGTCAGTTTGTGCCCTGCT
GTCCCTCCAGGACACTGAAATGCAGCGGATCTGCTCCTTACCACGTTTTTGAGGAATAAGAATCTTGAC
AGCTTCTCATGGCGCTGTTGACTATCTGTCTTATTACCTGCAGAGACTCTCGATGGAGAAGAATCCCC
AAAGCCGCGTGGTAGGCTGATAGAGAAGAAAGAGGTAGAGCTGGTTGTGAACAAGTTAGAAGATGCACA
GAAGTACCTGGCTCAGAAGTACTGCATCCTCGTCTGGGCTTGGGCATGGCGGATAAGCATCACTTGAGC
TGTGGAAGGAGAAAAATATCAGACACCCAAAAAGACTGGAATTTCTTGAGTCCTTTTATACCTTCTGTA
CCTGTATAGCCTGGATCGTCTTCCGACGTCACTACTTGACAGAGATCGAGGAAGAGGTAGGCAGGCTTTT
TCGGACCAATATGTTTAACATTCCTCGAAGGAAGCATGAAGACGAAGCGTCAGGAGGAGAGAAGAAGCGC
ATGACATTGGTACAGTTCAGGAGAATGATGGCAAAACGTCAGCAATTAAGGCGCGTGGACATGCGCT
CCCCAGTTTGTCTACGCTGCTGCCGTCTCTCAGAGAAAAAGCTCAACATATCACAGAGAAAAAGTATGT
GGCAGGCATCAAACTCCAACCCGAGAGGAGAAACATAATAACAGATCTGGAATCTGTGGCTATGCCAATA
GTTGGCATCTTGGGGGAGCCTCGAAATCTGTTCAACCCGAATACTCTCCTCCCCTGGAATTAGAAGAAA
ACACAAACCAAGTGGAAGAAGTTCTTCATAGTCGAGACAAATAACACAAGATCCAAACACACTGAG
TTTAGTCATGTCAAACTAGCTTCGCCAAGCACCTGCCCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>MR214302 representing NM_001163103 Red=Cloning site Green=Tags(s) MVKSEAMFTVPEFYSRKQFVGQSSTRLDQCGLRLGMWYWKDETKTLEFKSFTPAIELKEKGKKGKAVHF AEMDSGASERLTDKRFASRDEKAAKALEKRSQQGNVTLDDVKFVALLSLQDTEMQRICSFTTFLRNKNLD SFLMALLYYLSYYLQRLSMEKNPQSRVVGLIEKKEVELVVKLEDAQKYLAQKYCILVLGLGMADKHHLS CGKEKISDTQKDWKFFESFYTFCTCIAWIVFRRQYLTEIEEEVGRLFRTNMFMNIPRRKHEDEASGGEKKR MTLVQFRRMMAKRPAIKKAVDMRSPVLSTLLPSLREKAQHITEKKYVAGIKLQPREENIITDLESVAMPI VGILGEPRNLFNPNTLLPLELEENTKPSGRSSSIVETNNTKIQNTLSLVMSKLASPSTCPK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001163103
ORF Size:	1233 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001163103.1</u> , <u>NP_001156575.1</u>
RefSeq Size:	1453 bp
RefSeq ORF:	1236 bp
Locus ID:	210762
UniProt ID:	<u>D3Z0R2</u>
Cytogenetics:	12 C3
MW:	47.4 kDa