

Product datasheet for MR214325

Tmem52 (NM_027161) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem52 (NM_027161) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tmem52
Synonyms: 2310043I08Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR214325 representing NM_027161
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCCCGGGGCGTCAGCCACACAGGGGATCCTGCTGCTGCTGCCTCTCTGCCTCTGTCCCAGGTGA
 CGCTGGGTTCCGCGG[AT]CCGCAACTGTGACCCCTCGGATCAGTGCCGCTGTCCCCGAGGTGCCCGCC
 CCAGGCCCGCTGGAGCAGCCTATGGCATGTGGGGCTCATCCTGCTCGCCATACTCCTGATGCTTTTGTGT
 GGGGTACAGCCAGCTGTGTACGATTCTGCTGCCTCAGGAAGCAGACGCACACCCAGTCACACACGCCAG
 CGGCGTGGCAGCCCTGTGACGGGACAGTCATCCAGTGGACAGTGATAGCCCCGACACAGCACTGTGAC
 CTCTACAGTTCCGTGCAGTACCCACTGGGCATGCGGCTGCCCTGTACTTTGGAGAGCCAGACCCCTGAC
 TCCATGGTCCCTCCACCTACAGCCTCTATGCGTCTGAGCTGCCACCCTCCTACGATGAAGTTGTGAAGA
 TGATAAAAGCCAGGGAGGAAGTGGCAGCTCCCTCTGAGAAAACCAACTCTCTGCCTGAGGCCTTGGAGCC
 AGAGACCACTGGAGGGCCCCAGGAGCCAGGCCCCAGTCCCCAGCGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR214325 representing NM_027161
 Red=Cloning site Green=Tags(s)

MAPGPSATQGILLLLPLLPLSQVTLSAXXRNCPSDQCRCPRRCPPQARWSSLWHVGLILLAILLMLLC
 GVTASCVRFCCLRKQHTQSHTPAAWQPCDGTIPVDS SPAHSTVTSYSSVQYPLGMRLPLYFGEPPDP
 SMVPPTYSLYASELPSPSYDEVVKMIKAREEVAAPSEKTNSLPEALEPETTGGPQEPGPSAQR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/mm9025_c06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_027161

ORF Size: 609 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_027161.1](#), [NM_027161.2](#), [NP_081437.1](#)

RefSeq Size: 968 bp

RefSeq ORF: 609 bp
 Locus ID: 69671
 Cytogenetics: 4 E2
 MW: 22.3 kDa

Product images:



Circular map for MR214325