

Product datasheet for **MR214341**

Trif1 (NM_028081) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Trif1 (NM_028081) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Trif1
Synonyms: 2010012G17Rik; 4933421E11Rik; A1450568
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR214341 representing NM_028081
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCATCAATAGTAAAAAGGAAATTCAGAAAAAGGAAATAACAAGAAATATTCTCAAGGAAGCCATA
CTAATATTAAGGCATCATGTCTGAAGAATGATGCTGAATTTAAAAAGCTATTTGGTCTACTAAAGATT
GAGAGTATGTCTTACTCGAATTCCTGACCATCTGAGCTCTGAAAAAGTTTTAATTCCTTAACAGTTT
ATGAAAAGTAGTTCTTACAAAGATGCAACATTGGTGAAGAAAGAAAAAACAGAGTTTCTCTA
AGAAAAGAAAAGCGGAAACCATGAAGATGGGTAAACAAAAAGATAAAAAATTGAGAATGCTGATGACAC
GGTCATGAGCATCATGAATGGAAGTGTGCTAGTTCCCAACCTCTCAGTAGTATTTTACCAACATCA
GACATATCACAACATAACATTGTCAAGCCATAGTACAACCAGGGAAGACAAGAGAAGTGGCAGAAC
ACTGTTCTCATGAGAAGCAGGAAAAAGGCATTAAGTTCAAGCACATCTTTTGAACAAAGTACTTTTT
AAATAAAACTTTATGGAAGATATTTTCCAGTGACACCAGAACTAGAAGAAACCATCAGAGATGAA
AAAATAAGAAGACTTAAGCAAATTTAAGAGAAAAAGAAGCAGCTCTTGAAGAATTGCGTAAGAAGATGT
ACCAAAAACAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR214341 representing NM_028081
Red=Cloning site Green=Tags(s)

MASIVKKEIQEKGNKKYSQGSHTNIKASCLKNDAEFKKLFGLTKDLRVCLTRIPDHLSSGKSFNSFNLSL
 MKSSSYKDANIVVKKEKKQSFSKRKAETMKGNTKKIKIENADDTVMSIMNGTDVASSQPLSSILPTS
 DISQHNIVTSHSTTREDKRTEAEHCSHEKQEKGLSSSTSFEQSTFLNKNFMEDIFPVTPELEETIRDE
 KIRRLKQILREKEAALEELRKKMYQKQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_028081

ORF Size: 711 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028081.2](#), [NP_082357.2](#)

RefSeq Size: 1890 bp

RefSeq ORF: 714 bp

Locus ID: 321000


UniProt ID: [Q8CDD9](#)

Cytogenetics: 3 F2.3

MW: 27.5 kDa

Gene Summary: Together with SMCHD1, involved in chromosome X inactivation in females by promoting the compaction of heterochromatin. Also able to repress the ligand-induced transcriptional activity of retinoic acid receptor alpha (RARA), possibly through direct recruitment of histone deacetylases.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR214341