

Product datasheet for MR214395L4

Pcdhb5 (NM_053130) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pcdhb5 (NM_053130) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Pcdhb5
Synonyms:	Pcdhb4A; PcdhbE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214395).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_053130
ORF Size:	2376 bp



[View online »](#)

