

Product datasheet for **MR214426L3V**

Swsap1 (NM_025870) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Swsap1 (NM_025870) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Swsap1
Synonyms:	2310047B19Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025870
ORF Size:	834 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214426).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025870.1 , NP_080146.1
RefSeq Size:	1846 bp
RefSeq ORF:	837 bp
Locus ID:	66962
UniProt ID:	Q8VCI7
Cytogenetics:	9 A3
Gene Summary:	ATPase which is preferentially stimulated by single-stranded DNA and is involved in homologous recombination repair (HRR). Has a DNA-binding activity which is independent of its ATPase activity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)