

Product datasheet for **MR214440A1V**

Mouse Olfr1306 (NM_001011803) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: Olfr1306

Synonyms: MOR245-15

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR214440 representing NM_001011803
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAAGAGTGAATCAGTCTGTGGTGTCTCAGAGTTTGTGTTCTGGGACTCACCAACTCCTGGAGTATTC
AACTATTTCTCTTTGTGTTCTCCTCCACGTTTTATGTGGCAAGCATGATGGGAAACAGCCTCATCGTGTT
CACTGTGGCATCTGACCCTCACTTACACTCTCCATGTACTTTCTTTGGCTAACCTCTCTTTCATTGAC
TTGGGTGTTTCTTCTGTTACTTCTCCAAAGATGATTTATGACTTATTAGAAAGCACAAAGTCATCTCCT
TTACAGGTTGTGTTACTCAAATATTTTCATTCATGTCATTGGTGGTGTGGAGATGGTGCTGCTCATAGC
CATGGCCTTTGACAGATATGTGGCCATATGTAAGCCTCTCCATTACCTGACCATTATGAGCCCGAGGATG
TGCATCTTGTTTTAGTGGCTGCCTGGGTGGTGGCCTTATGCATTCTTTGATTCAACTGGCTTTTTTAG
TAAACTTACCATTCTGTGGACCAAATGTGTTGGACAGCTTCTACTGTGACCTTCCTCGGTTTATCAAAC
AGCCTGCATAAATACTTACCAGCTGGAATTCATGGTCACAGCCAACAGTGGGTTCAATTTCTGTGGCTCC
TTCTTCATACTGATCATTTTATATATTGTCATCATAATCACTGTGCAAAAACATTCTCAAGTGGTTCCT
CTAAAGCCCTATCTACTCTTTCTGCTCACATCTCTGTGGTGTCTTATTCTTTGGTCTTTGATATTTT
CTATACCTGGCCATCTCCTTCAACACATCTTGATAAGTTTCTGGCCATATTTGATGCAGTTATCACTCCG
TTTCTGAATCCTATAATCTATACATTTAGAAATCAAGAAATGAAGGTAGCAATGAAGAGAGTATGCAGAC
AGCTAGTTAGCTATAGGAAAATTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>MR214440 representing NM_001011803 Red=Cloning site Green=Tags(s) MERNQSVVSEFVFLGLTNSWSIQLFLFVFSSTFYVASMMGNSLIVFTVASDPHLHSPMYFLLANLSFID LGVSSVTSPKMIYDLFRKHKVISFTGCVTQIFFIHVIGGVEMVLLIAMAFDRYVAICKPLHYLTIMSPRM CILFSVAWVVGMLMHSIQLAFLVNLPCGPNVLDSEYCDLPRFIKLACINTYQLEFMVTANSQFISVGS FFILIIISYIVIIITVQKHSSSGSSKALSTLSAHISVVFLFFGPLIFFYTWSPSTHLDKFLAIFDAVITP FLNPIIYTFRNQEMKVAMKRVCRQLVSYRKTS TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001011803
ORF Size:	936 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001011803.2, NP_001011803.2</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	258023
Cytogenetics:	2 E3
MW:	35.2 kDa