

## Product datasheet for **MR214487**

### **Olfr1226 (NM\_146967) Mouse Tagged ORF Clone**

#### Product data:

|                           |                                                                             |
|---------------------------|-----------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                         |
| Product Name:             | Olfr1226 (NM_146967) Mouse Tagged ORF Clone                                 |
| Tag:                      | Myc-DDK                                                                     |
| Symbol:                   | Olfr1226                                                                    |
| Synonyms:                 | MOR233-2                                                                    |
| Mammalian Cell Selection: | Neomycin                                                                    |
| Vector:                   | pCMV6-Entry (PS100001)                                                      |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                        |
| ORF Nucleotide Sequence:  | >MR214487 representing NM_146967<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTTAAACCACAGTAGTGCTACTGAGTTCATTCTTCTTGGACTTTCACAGAACTCAAAGTTGAGAAAG  
 TACTGTTTGTATATTTTGTGATCTACCTTGCAACTATTGGGGCAACATGATTATTGTAGTAACCAT  
 CATCTATAGCCCTGCACTGCTGGGCTCCCAATGACTTCTTCTTAATATTCCTGTCCTTCTGGATGCT  
 TGCACTTCTTCCACTGTACCCCAAGATGATTGTAGACTTCTTCTATGAAAAGAAGACCATCTCCTTTG  
 AATGCTGCATGACAACTGTTTGTGCTGCACTTCTTCACTGGGATGGAGGTGATTGTCCTTTCAGCCAT  
 GGCTATGACCGCTATGTGGCCATTGCAAGCCTCTTCACTACTCCTCCATCATGACCAGGAGACTCTGT  
 GGCATTTTGGTGATGTTGTCTTGGGCAGGAGGCTTCTTGCACTTCTATCATACAAATCATCTTCACTTTGC  
 AGCTGCCCTTCTGTGGACCAATGTCAATTGATCATTACATGTGTGACTTGTCCCACTACTGAAGCTTGC  
 CTGCACTGACACACATATTTGTCAATTTAGTGTGTTGCCAAGTGGGTCTATCTGCATTATAATCTTC  
 TCTATTTTGTCTGCTCCTATAGTGTCTTGTCTCTCTGAGAGCCCACAGTTCTGAAGGGCGGCGTA  
 AAGCTCTCTCCACCTGTGGATCCACATTACAGTTGTGCTTTTGTGTTTGTGTCATGCATTAATATA  
 TGCTCGACCTGCTGCCTTCTTCTGAAAAAATGCACTTGATTTGCCACAATCATAACACCACCTCTG  
 AATCCTATGGTATATACTTTCAGGAATAGGGAATGAAGAATGCTATTAGGAACTGTGGAAGAGATGGA  
 AAGTGTTTCTGGTGACATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

**Protein Sequence:** >MR214487 representing NM\_146967  
 Red=Cloning site Green=Tags(s)


MLNHSSVTEFILLGLSQNSKVEKVLFIYLLIYLATIGGNMIIIVVTIIYSPALLGSPMYFFLIFLSFLDA  
 CTSSTVTPKMIVDFFYEKKTISFECCMTQLFAVHFFTGMIEVIVLSAMAYDRYVAICKPLHYSSIMTRRLC  
 GILVMSWAGGFLHSIIQIIFTLQLPFCGPNVIDHYMCDLFPLLKLACTDTHIFVILVFANSGSICIIIF  
 SILLVSYSVILFSLRAHSSEGRRKALSTCGSHITVVLLFFVSCILYARPAAFSSEKNALVFATIIITPLL  
 NPMVYTFRNREMKNNAIRKLWKRWKVVSVDI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9069\\_g10.zip](https://cdn.origene.com/chromatograms/mm9069_g10.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_146967

**ORF Size:** 930 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_146967.1](#), [NP\\_667178.1](#)

**RefSeq Size:** 933 bp

**RefSeq ORF:** 933 bp


**Locus ID:** 258969

**Cytogenetics:** 2 E1

**MW:** 35 kDa

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

## Product images:



Circular map for MR214487