

Product datasheet for **MR214510L3V**

Fam83b (NM_001045518) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fam83b (NM_001045518) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fam83b
Synonyms:	C530008M07Rik; Gm516
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001045518
ORF Size:	3036 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214510).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001045518.1
RefSeq Size:	3136 bp
RefSeq ORF:	3039 bp
Locus ID:	208994
UniProt ID:	Q0VBM2
Cytogenetics:	9 D
Gene Summary:	Probable proto-oncogene that functions in the epidermal growth factor receptor/EGFR signaling pathway. May activate both the EGFR itself and downstream RAS/MAPK and PI3K/AKT/TOR signaling cascades.[UniProtKB/Swiss-Prot Function]



[View online »](#)