

Product datasheet for **MR214543A1V**

Mouse Olfr895 (NM_146875) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr895 (NM_146875) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr895
Synonyms:	MOR170-1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR214543 representing NM_146875 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCACCTCAATGCACAAAAGACTATGGAAATGATTCTTCTGTGTCTGAGTTTATTCTTATGGGACTGA
CAGACCAACCTGAGCTCCAGCTGCCCTTATTTGTTCTGTTCTTGGTGAACACAGTCACTGTGATGGG
AAATTTGAGCTTAATGAATCTCATTGCTTAACTCAAACCTTACACCCCATGTACTTTTTCATCTTC
AACCTGTCCTTCATTGACTTCTGTTATTCAATGGTCTTTACCCCAAAATGCTGATGGGCTTTGTTGTAG
AGAAGAACATCATCTCCTTCAGAGGATGCATGACTCAGCTGTTTTCTTCTTATTTTTGTGAACCTCGA
GAGTTATGTGCTGACAGCAATGGCCTATGATCGCTATGTGGCCATCTGTGAGCCACTGCTCTACAAAGCT
GTTATGTCACCTGGCATCTGCTTTCTATTGATATTTGTACTTATTTGATGGGGCTTGTTAGTGCCTTAT
TCCATACAGGTTTTATGATAAGGCTCAACTTTTGTGATTCTAATGTTATCAACCACTATATGTGTGATAT
CTTCCCTCTCTTCGGCTCTCCTGCAGTAGCACCTATCTCACTGAGCTTGTGAGCTCTGCTGTGGTGGGT
ACAGCTATTATCTTATGTTGCCTCATTATCTTAATCTCTTATGGTATGATCCTTTACAATATCATTATA
TGTCATCAGGTAAGGGTTGGTCCAAAGCCTTGGGCACTTGTGGTCTCACATTATACTGTTAGCTTGT
CTATGTAACCTGGGATGCTTGCTTATGTCAAGCCATCATCTGCTGAGACTGTGGGACAGGGGAAAAATTTTC
TCAGTATTTTACACTTTTTTGGTGCCCATGCTGAATCCTCTTATTACAGCCTCAGAAACAAGGATGTCA
AACTTGCTGTGAAAAAACCTGGAAGAGACTCATATGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR214543 representing NM_146875 Red=Cloning site Green=Tags(s) MSLNAQKTMENDSSVSEFILMGLTDQPELQLPLFVLFVNYTVTVMGNLSLMNLICLNSNLHTPMYFFIF NLSFIDFCYSMVFTPKMLMGFVVEKNIIISFRGCMTQLFFLFFVNSESYVLTAMAYDRYVAICQPLLYKA VMSPGICFLLIFCTYLMGLVSALFHTGFMIRLNFCDSNVINHVMCDIFPLFRLSCSSTYLTELVSAAVVG TAIILCCLIIILISYGMILYNIHMSGKGWSKALGTCGSHIITVSLFYVTGMLAYVKPSSAETVGQGKIF SVFYTFLLVPMNLPLIYSLRNKDVKLAVKKTWKRLTC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146875
ORF Size:	948 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146875.2</u> , <u>NP_667086.2</u>
RefSeq Size:	951 bp
RefSeq ORF:	951 bp
Locus ID:	258875
Cytogenetics:	9 A5.1
MW:	35.7 kDa