

Product datasheet for **MR214610**

C2cd2 (NM_174847) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C2cd2 (NM_174847) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C2cd2
Synonyms:	5730563M15; 5830404H04Rik; AI314146; AI449533; BF168443; mKIAA4045; ORF25
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR214610 representing NM_174847
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGGCTCGTGGCTCAGCGAGGTGCAGTGGCTCTTCTGGTGTGCTCTTCGTCGCGGCCCTGGGCACCG
TGGCCTGTACTTGGCGCAGTGGGCGCTGGCCAAGGCGCGACCCCGCCCGGCGGCGGGGAGCCGGA
TGAGCTAAGGCGCCGGGAGTCCGACACACTCCTTTCTGGATACTGACACGGGATAGCTGGGGAAACCG
TGGCAGCGGCTTGGGTGACCGCCTGAACTATGAGGCCGAAAAGAGAGGGGGTCCACTGCGTCTCTCTT
TCCAGAAGGACCCACGGCCACAATCCTTGACGCTGACAGTTGAGAAGGTATCCAGTGTAGTCAAGTCAAC
CCAGGAGAAGGTGGTGTCTGTCTGTGGTGGGTGAGACCCTTCAGTTCTTGGTCACTGACGGCCCTGCC
TCTGCCACAGGCTCGGAGTGCAGCTGTATGATGCCATCTTTCCCATTTTCAATTAAGGTAGAATTCC
ACATGGAAGAGAAGAGAGAGGACATACAAATCAGGTGGTCTTTCACCCACGTCCCAGAAACGCCATCAA
GATCCAGCCCCAGGCCAGGGGAGAAGCAGGCGCTCGGGGTAACATGCTTTCTGAAGCCCTTGAGGAT
CTTTTCAAGCACCTAGTTAATGCTGCCTCTCCATCTGTTTTCTGAGCACCAAACCCACGAGGTCAAAG
AAGCTCAGAGTCTACAGTGCCCTTCTCCACTGCTCAGGAACCCCTGCCCTCCCAAACCTCCAAGGGCTCA
TGAATAAAATTCAGGTGAAGAATATCAGAGTTTCCCTGATAAACCCACCTGGTCTTCAGGCCTTAGC
CATGTGTGCGTGGCTCAGCTGAACGATCCTGAACAGAGATTATCAGCACCCCTAGTGAGAAACACTACTG
ACCTCAGTTGGGAGGAGGAATTCACCTTGAAGTGAACGCCAAATCCAAGGAGCTGGTCTGCAGATCTC
ACAGGATGGGTGCTCCTCAGAAGGTCTGTTGGGGATAGCAACCATTACCTGGACCTATTTAGGAAGCAG
CCTAATGGCCACAGACCTTCAGACTGATCAGTGGGACAGAGCCTGACAGCTTGGTGTGGGCTCAGTCA
CTGCAGAGTTTTCTTACGTAGAACCTGGGAGCTGAAATCCTGGCCAGCTCCTCCCCTGTTTTCTGCTGC
AAAGATAGAAAAGGACCGTACGGTGATGCCATGTGGGACTGTGGTCAACCTGTCACTGCTGTGAAGACC
AAACCTCGGTTTGACACTGGGAGGGCAACCCACTGAACTCAGAGTCCCCTGTGAGGACACCCGTCACAG
TGAAGGTCATTGAGAAGGACATCTCTGTCCAGGCCATTTCTTGCCACAGCGCCCCGTCAGCAAGACATT
CTCCTCTCAGACACGGAAGTGTGGTGTGAATGGCTCAGACCCTGTGGCTGAAGTGCATTCGGCAA
CTCAGCGAGTCTCAAAGCTGAACTCAAGTCCCACGGAAGAAAAGCACCATCATAATTCAGGGATCT
CCAAGACCTCACTATCTCAGGACCACAATGCTGCCCTCATGCTGGACTACGAGCTTCTATGGACAGCAC
CAACCAGGAGATGCCACGTCTGCTTTGTGCCATCCAGAGGCCACTGAAGCCTCTGCTACTACCCACCT
GAAGAAAATGAGCCAGCCAGACCCTACCTGCCCTCAAGCCTAGAGAAAATGATCTAGATTCTGGGAAC
TGGAGAAGGAGTCTCCAGTAGCATCATGGAGCGGACCAGCCTGCAGGAGCCTGACGGTGTGAAGTATC
CGAGAGTCCCTGAGTACTTCAGAGCTAGGAGCTATGAAGAAGCATAAAGGAGGACTCCTAAGGAAAGGT
GCCAAGCTCTTCTCCGCCGGGACACCAACAGAAAGACCCGGGCTTGAGCCAGTACATAATGACCTGG
TATTCCTGCAGCAGCCAGAGGGGAGGCAGAAGAAGGGGGCAGCCCTCGGCCGGCTCTGAACAGAAAGCT
GCTTACCAGGCACAGAGGCAAACACACTATGAATGGAGTACCCAGAGAGCCCTGCATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR214610 representing NM_174847
Red=Cloning site Green=Tags(s)


MGSWLSEVQWFLVLSLFVAALGTVGLYLAQWALAKARPPRRRAEPDELRRRESDTLLSWILTRDSWGNQ
 WQAAWVTALNVEAEKRGGLRLSFQKDRPQSLQLTVEKVSSVVKSTQEKVVICHVVGETLQFLVSAGPA
 SATGSECQLYDVHLSPFHLKVEFHMEEKREDIQRWSFTHVPETAIKIQPQAPGEKQALGVNMLSEALED
 LFKHLVNAASPSVFLSTKPTQVKEAQLQCPSSSTAQEPCKPPRAHELKLQVKNIRVSLINHPGASGLS
 HVCVAQLNDPEQRFISTLVRNTTDLSEWEEFTFELNAKSKELVLQISQDGCSEGLLGIATIHLDLFRKQ
 PNGPQTFRLISGTEPDSLVLGSVTAEFYVEPGELKSWPAPPVSAKIEKDRTVMPCGTVVTTVTAVKT
 KPRFDTGRATPLNSESPVRTPVTVKVEKDISVQAISCHSAPVSKTFSSDTELLVLNGSDPVAEVAIRQ
 LSESSKLLKSPRKSTIIISGISKTSLSQDHNAALMLDYAASMDSTNQG DATSALCHPEATEASATTPP
 EENEP AQTLPALKPRENDLDSWELEKESPVASWSGPALQEPDGDDELSESSLSTSELGAMKKHKGGLLRKG
 AKLFFRRRHQQKDPGLSQSHNDLVFLQQPEGRQKKGATLGRLLNRKLLTRHRGKHTMNGVPREPCI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9013_g09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_174847

ORF Size: 2088 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_174847.2, NP_777272.2</u>
RefSeq Size:	6568 bp
RefSeq ORF:	2091 bp
Locus ID:	207781
UniProt ID:	<u>E9Q3C1</u>
Cytogenetics:	16 C4
MW:	77 kDa

Product images:


Circular map for MR214610