

Product datasheet for **MR214634**

Olfr362 (NM_147051) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr362 (NM_147051) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr362
Synonyms:	MOR158-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR214634 representing NM_147051 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGAGCTGTGCCCGAATGCTTCACATTCTCCAATCTTCTTGCTCCTGGGTCTCAAGGGCTGGTGTGC
 CCCACACTTTCTCTTTTATTGTTCTATTTATTTATCTGACCACCATATTGGGAATGTGACTTTGGT
 ACTGCTCATCTCTGGGACTCCAGGCTCCACTCACCTATGTATTATCTCTTGGTGGCTCTCCATGATA
 GACTTGGGACTATCCACAGTCACCCTGCCACAGCTGTTGGTTCATCTGTCATCTGATTCTCCAGCCATTC
 CTGCTGCCCGCTGTTTGACTCAGTCTTTTTCTTTATGCATTTGGAGTCACAGATACACTTGTCTATTGC
 TGTCTATGGCTCTGGATCGCTATGTGGCCATCTGTGACCCCTCTGCACTATGCTTTGGTGATGAATCGACAA
 ATTTGTGCCCGGTTACTGGCCTTGAGCTGGGTGGTATCTATAGTACACACCATGTGTCATGTGGGACTCA
 TCCTGCCACTGTGTTGGGCTGGAGATGCTGGGGCAATGTTAACTTCCTCACTTCTTTTGATACCCG
 ACCACTCCTGCGAGCCTCTTGCTCTGACACACATTCCAATGAGCTGGCTATATTCTGGAGGGTGGCTTT
 CTAATGTTGGGACCTGTTCCCTTATTGTACTTTCCTATGCCGAATTGGGATCACCATCTTACGATTGC
 CTTGAGCTGCTGGTCGCCGACAGCAGTCTCCACCTGTGGATCCACCTCACCATGGTTGGCTTCTGTGA
 TGGCACAATCATTTGGGTCTATTTCCAGCCTCCCTCACAGAACTCCCGGAATCAAGACATGGTGGCTTCA
 GTAATGTACACTGCTATTACCCATTGGCTAACCCATTTGTTTATAGTCTCCGTAATAAGGATGTGAAAG
 GTGCCCTCCACAGACTACTGAGACAGGGGAGAGTGGACTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR214634 representing NM_147051
 Red=Cloning site Green=Tags(s)

MSCAPNASHSPIFLLLGFSRAGVPHTFLFLFLFIYLTITLGNVTLVLLISWDSRLHSPMYLLRGLSMI
 DLGLSTVTLPLQLLVHSSDSPAIPAARCLTQFFFFYAFGVTDTLVIAMALDRYVAICDPLHYALVMNRQ
 ICARLLALSWVVSIVHTMLHVGLILPLCWAGDAGGNVLPFFCDHRPLLRASCSDTHSNELAIIFLEGGF
 LMLGPCSLIVLSYARIGITILRLPSAAGRRRAVSTCGSHLTMVGFLYGTIIWVYFQPPSQNSRNQDMVAS
 VMYTAITPLANPFVYSLRNKDKGALHRLLRQGRVDS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9021_d11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_147051

ORF Size: 951 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_147051.1](#), [NP_667262.1](#)

RefSeq Size: 954 bp

RefSeq ORF: 954 bp


Locus ID: 259053

Cytogenetics: 2 B

MW: 35.5 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR214634