

Product datasheet for **MR214719L4V**

Ap5b1 (NM_001033448) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ap5b1 (NM_001033448) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ap5b1
Synonyms:	Gm962
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001033448
ORF Size:	2628 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214719).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033448.2 , NP_001028620.1
RefSeq Size:	3092 bp
RefSeq ORF:	2631 bp
Locus ID:	381201
UniProt ID:	Q3TAP4
Cytogenetics:	19 A
Gene Summary:	As part of AP-5, a probable fifth adaptor protein complex it may be involved in endosomal transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)