

Product datasheet for **MR214730**

4921501E09Rik (NM_001009544) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	4921501E09Rik
Synonyms:	PHF8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >MR214730 representing NM_001009544
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCCTTGGTGCCTGTGATTGTCTTTGTCGACAGCCTTATAATGTGAACCACTTCATGATTGAGTGCG
 GCCTGTGCCAGGACTGGTTTCATGGTAGCTGTGTTGGGATTGAAGAAGAAAACGCTGTTGACATTGACAT
 CTACCACTGCCAGACTGTGAGGCTGTATTTGGCCTTCCATTATGAAAACTGGCATGGATCGCCTAGA
 GGGCATGATGCTCTCAAGAGGAAGCCATTGATGACTGGGAGCCCTACCTCTATTCACCAAATCCAGGGTA
 GGATTTTGAAGCTCTGACAAAGTGATTTGAAGCCACAGGAAGTCATTTAACAGTAGAATTTCTGCT
 AGAAAGTAGCTTCAGCATACCTATCCTTGTCTTGAAGAAGGATGGGTTGGGGATGACACTTCCACCGACA
 TCATTCAGTGTGAGAGATGTTGAGCAATATGTTGGTCTGACAGAGAGGTTGATGTGATTGATGAGCCT
 GTCAGGCTAGCAGCAAGATGATGTTGGGAGATTTGTGAAATATTATTACAGCAGGAAGAGGGAAAAGAC
 CTTTAATGTTATTAATTTGGAATCTCTAATACTGGGCTTTCTAATATAGTAGAGATGCCCAAGATAGTT
 CGTAAGCTTTCTGGGTTGAAAACCTATAGCCAGAAAGATTGCGGCTTTGAGAGACCCAATGTACAGAAAT
 ATTGCTTATGAGTATGCAGAATAGCTATACAGATTTTACATTGAATTTGGTGGCACTTCAGCCTGGTA
 CCATGTAGTTAAGGGTGAGAAGATCTTCTATCTGATCCGCCAACAGATGCTAACCTGACTCTTTTGTAG
 TGTGGAGCAACTCTTCAACCAAGGAAATATTCTTTGGTGACCAGGTAGACAGATGTTACAAGTGTT
 CTGTGAAGCAAGGACAGACACTTTTCACTCTACAGGATGGATCTATGCTGTGCTGACCCAGTGGATTG
 CCTAGCATTGGGGGAACTCTTACACAGTTTAAACATTGAAATGCAGCTCAAAGTCTATGAGATTGAG
 AAGCATCTGTGTACAACAGACCTCTCAAGTTTCAAGCTTTGAGACCATCTGTTGGCATGTGGGAAAC
 ACATCTGGACATCTTTAAGGTTTACGGGAGAATGGGAGATGCTCTGCATCTTCTTGGTCCATGGTGG
 TAAAGCTCTAACTTGGCCTTTCAAGCCTGGACAAGGAAAGAGTTCTGCCTAACCGTGATATGAGATC
 CCAAAGACAGTGCAGACCATACAGCTCATTGAAGATCTGGCTTGGGAAATTCATCTGGTTGAAGATAATT
 TCCAACAGAACATCCCAAGACAAGTTCTACCTTTAAGTTCAGCAGCTCCATACCAATTCCTCTTCCTT
 AACTAGGCAAGCTCATCCCACTTCAGTATCCAAGTCTAAATATCACTGCCCTCCAAAAATGGTTCAGAA
 AAGAAAAGCCCAAAGTATAAGACACCTTTAAGAAGGAAGAACAATGAGCAAGGAGAGTTCAAGCTCTGG
 GGCCTAATTATCAACATAGCTGTAATTTCAAGGATCTAAATAATGATCAACAATGGAAGGTAGGCTGTCT
 CCAAAACGTAGAGTTTAACATCACTAGTACTTGTGTAATGATTCAGATGATAACTCAAAACATGACAAA
 GCCTTCAATGGGAATGAGAGCCTGGTGCCTCTGTTGATGGCTAATGGCAATACGAAGAGGGAAAAAGCT
 TGTCCGAATCTTGGCAAGCCAAAATCACAAGAAGGAAGATACATCAAGGTTGATGAAAGAACAGCGAT
 GGTAGACACATTTGACTTGGATTTAAATGATGAGCTACAGATTGAAAAGGGATTAGGTAAAGAGAAAGTG
 AGTCTGGTAGTAAGACCAAAAATTCCTCCAAGTTGCCCATGTGAAGCTTTGCTCTGACCAAAACCAAA
 GTCATGATCCAGGGGAACTGAGATCACCATTGAGAGAACTATAAAATAGATAAGGAAATGGTGAAGA
 GGTGGAAGGCAAGTTTAAGAATGGAAGTGAATCTGGTGGAACTATTGATCCACTAAAGGCCAAGAAACAG
 GTGGGTGGACCTGGCTGTGTTATCCTCACTGATGACCCAACTCTCCAGCAATAAGAAGGACATTCAGG
 ACATGTGTGCATGTCTGATTGTCAGTCTTCTCCTCCTTCAAGCTCCCTCCAGCCTTCAGGCATTTTG
 GAGTGGGGACAAGACCAAGGCAACAAAAGCCCTAGCAGTGTCTTGGGCACAGTGTCTGACAGTCTGTT
 TCTCGGCGTACCCAAGGAAATGGCCCATCAAGCGGCTAGCATACTGGAGAAATGAAAGTGAGGCGGGAA
 ATGCTTGTGTAGATGAACAGGGCAGCTTAGGAGCATGCTTCAAGGAAGCTAAGTATATCTTCTTATCCCT
 GGAATCTGATGACGACGACGATGATGATGATGATGACGACGATGACAACAACCTGCTTTGAGAGTTCTA
 CCAAGAAAAAGGAGATATCAGATAATGCCCAAGGATTCTAGGGCAGATGAGATCCCAACACTACCAA
 AGCAGGACCATTTCTGTGCTGAGGGTAACAGAGTTGGCTCCATTGAGATAGACTGGACTACAGCAGCTAC
 AAAGCTGGCACAGGAGGTGAACAACATGAAGCTCACTCTCACTAGTCCAGGAAAAAGCAGCCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR214730 representing NM_001009544

Red=Cloning site Green=Tags(s)

MALVPVYCLCRQPVNVNHFMIIEGLCQDWFHSCVGIIEENAVDIDIYHCPDCEAVFGPSIMKNWHGSPR
GHDALKRKPLMTGSPTSIIHQIQGRIFESSDKVILKPTGSHLTVEFLLESSFSIPILVLKKGDLGMLTPPT
SFTVRDVEQYVGSREVDVIDVACQASSKMLGDFVKKYYSRKREKTFNVINLEFNTGLSNIVEMPKIV
RKL SWENLWPEDCGFERPNVQKYCLMSQNSYDFHIEFGGTSAWYHVKGKIFYLIRPTDANLTLFE
CWSNSSNQKEIFFGDQVDRYKCSVKQGQTLFIPTGWIYAVLTPVDCLAFGGNFLHSLNIEMQLKVYEIE
KHLCTDLFKFSPFETICWHVGHILDFKGLRENGRCPASFLVHGKALNLAFAWTRKEVLPNRDNEI
PKTVQTIQLIEDLAWIEHLVEDNFQNIPTSSTFKFQQLHTNSSSLTRQAHTSVSKSLSLPSKNGSE
KKSPKYKDTFKKEEQMSKESALGPYQHSCNFMDLNNDQWQVGLQNVFNITSTCLNDSDDNSKHDK
AFNGNESLVPLL MANGNTKREKSLSESWQAKITKKEDTSRLMKEQAMVDTFDLNDELQIEKGLGKEKV
SLVVRPKIPPSLPHVKLCSDQNQSHDPGETEITIVENYKIDKEMVEEVEGKFKNGSESGGTIDPLKAKKQ
VGGPGCVILTDDPNSSNKKDIQDTCMSDLQSSPPSQAPSSLQAFWSGGQDQGNKSPSSVLGTVDSPV
SRRTPRKWPIKRLAYWRNESEAGNACVDEQGSGLGACFKYIFLSLESDDDDDDDDDDDNNPALRVL
PKKRKISDNAPRIPEADIEPTLPKQDHSVPEGNRVGSIEIDWTTAATKLAQEVNMMKLTLTSPGKSSQ

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mm9096_d12.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:




ACCN:	NM_001009544
ORF Size:	2724 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001009544.2 , NP_001009544.1
RefSeq Size:	3916 bp
RefSeq ORF:	2727 bp
Locus ID:	74042
UniProt ID:	Q80TJ7
Cytogenetics:	17 B1
MW:	102.1 kDa

Gene Summary:

Histone lysine demethylase with selectivity for the di- and monomethyl states that plays a key role cell cycle progression, rDNA transcription and brain development. Demethylates mono- and dimethylated histone H3 'Lys-9' residue (H3K9Me1 and H3K9Me2), dimethylated H3 'Lys-27' (H3K27Me2) and monomethylated histone H4 'Lys-20' residue (H4K20Me1). Acts as a transcription activator as H3K9Me1, H3K9Me2, H3K27Me2 and H4K20Me1 are epigenetic repressive marks. Involved in cell cycle progression by being required to control G1-S transition. Acts as a coactivator of rDNA transcription, by activating polymerase I (pol I) mediated transcription of rRNA genes. Required for brain development, probably by regulating expression of neuron-specific genes. Has activity toward H4K20Me1 only when nucleosome is used as a substrate and when not histone octamer is used as substrate. May also have weak activity toward dimethylated H3 'Lys-36' (H3K36Me2), however, the relevance of this result remains unsure in vivo. Specifically binds trimethylated 'Lys-4' of histone H3 (H3K4me3), affecting histone demethylase specificity: has weak activity toward H3K9me2 in absence of H3K4me3, while it has high activity toward H3K9me2 when binding H3K4me3.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR214730