

Product datasheet for **MR214743L3V**

Tmem8b (NM_001085508) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem8b (NM_001085508) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem8b
Synonyms:	4930500O05Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001085508
ORF Size:	1416 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214743).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001085508.2 , NP_001078977.1
RefSeq Size:	4787 bp
RefSeq ORF:	1419 bp
Locus ID:	242409
UniProt ID:	B1AWJ5
Cytogenetics:	4 A5-B1
Gene Summary:	May function as a regulator of the EGFR pathway. Probable tumor suppressor which may function in cell growth, proliferation and adhesion (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)