

Product datasheet for **MR214744A1V**

Mouse Olfr952 (NM_146503) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Olfr952 (NM_146503) AAV Particle
Tag: Myc-DDK
Symbol: Olfr952
Synonyms: MOR171-21
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR214744 representing NM_146503
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAAGAATTAATCATACCCAGTGGCAGAGTTCATTCTTGCTGGTTAACAGAGAATCCAGAACTCC
 AATTGCCTCTATTCTTATCTTCCTATCTGTCTATTTGTTTACAGTTGTTGAAACTTAGGAATGATTGT
 CTTGATTCTGATTAGTTCTCAGCTACATACACCTATGTATTATTTACTCAGCAGTCTGTCCTTTATTGAC
 TGCTGTCAGTCTACTGTCATTGTTCTCTAAATGCTGTTGAACTTTGTGACAGAGAAGAATGTCATCCTGT
 ACCCAGAATGCATAGCTCAGTTCTATTTCTTCTGTACTTTTGTGTTGCAAAATGTCACATGTTGGCTGC
 AATGGCATATGACCGCTATGTGGCTATATCTAACCCCTTTGCTTTACAAAGTAACCATGTCCTATCAAGTC
 TGTTTGTGATGGTAGCTGTGGTGTATGGTATTGGCTTAATCAGTGCCACAGCTCACACAGTCTTCCTGC
 TAAGATTGTTTTCTGTAAGGCTGATAAAATAAACCACTACTTTTGTGATCTTTTCCATTACTTGAGCT
 GTCTTGCTCTAGTACTTTTATCAATGAAATATTAGCACTGTCCTTCAGTGCATTTAATATTATTGTACCA
 GCCATGACCATCATTGGCTCTTACATCTTCATCATTATCAGCATTTTGCACATTAATCCTCTGGGGCA
 GAGTCAAGGCCTTAGCACCTACAGCTCCCATATCTTGCTGTTGCTATTTTCTTGGTTCTACAACATT
 CATGTATCTACAACCATCTTCAGTCAGCTCCATGGACCAATGGAAAGTGCATCTGTGTTTATACCATT
 GTTGTGCCTGTGTTGAATCTCATGATCTACAGCCTGAGAAATAAAGATGTCAAAGTTGCCCTAAAAATGT
 TACTACAGAAGATGTTTCCTCAAAACAAGGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>MR214744 representing NM_146503 Red=Cloning site Green=Tags(s) MEELNHTPVAEFILAGLTENPELQLPLFLIFLSVYLFVVGNLGMIVLILISSQLHTPMYYLLSSLSFID CCQSTVIVPKMLLNFTVTEKNVILYPECIAQFYFFCTFVVAKCHMLAAMAYDRYVAISNPLLYKVTMSYQV CLLMVAVVYGIGLISATAHTVFLRLFFCKADKINHYFCDLFPLLELSCSSTFINEILALSFSAFNIIVP AMTIIGSYIFIIISILHIKSSGGRVKAFSTYSSHILAVAIFFGSTTFMYLQPSSVSSMDQWKVSSVFYTI VVPVLNLMYSLRNKDVKVALKMLLQKMFPQNKE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146503
ORF Size:	942 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146503.1, NP_666714.1</u>
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	235248
Cytogenetics:	9 A5.1
MW:	35.4 kDa