

Product datasheet for MR214746L3

Smpdl3a (NM_020561) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Smpdl3a (NM_020561) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Smpdl3a
Synonyms:	0610010C24Rik; AI529588; ASM3A; ASML3; ASML3A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214746).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_020561
ORF Size:	1335 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_020561.2 , NP_065586.3
RefSeq Size:	1791 bp
RefSeq ORF:	1338 bp
Locus ID:	57319
UniProt ID:	P70158
Cytogenetics:	10 B4
Gene Summary:	Has in vitro nucleotide phosphodiesterase activity with nucleoside triphosphates, such as ATP (PubMed:26792860). Has in vitro activity with p-nitrophenyl-TMP. Has lower activity with nucleoside diphosphates, and no activity with nucleoside monophosphates. Has in vitro activity with CDP-choline, giving rise to CMP and phosphocholine. Has in vitro activity with CDP-ethanolamine. Does not have sphingomyelin phosphodiesterase activity (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR214746L3