

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for MR214766L4V

Magea3 (NM_020017) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Lentiviral Particles
Magea3 (NM_020017) Mouse Tagged ORF Clone Lentiviral Particle
Magea3
Mage-a3
Puromycin
pLenti-C-mGFP-P2A-Puro (PS100093)
mGFP
NM_020017
963 bp
The ORF insert of this clone is exactly the same as(MR214766).
The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<u>NM 020017.2, NP 064401.2</u>
1107 bp
963 bp
17139
X 72.38 cM



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US