

Product datasheet for **MR214842A1V**

Mouse Olfr110 (NM_146328) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr110 (NM_146328) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr110
Synonyms:	MOR249-2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR214842 representing NM_146328 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAGGAAAGAATCAAACAGCTCCATCTGAATTCATCATCTTGGGGTTCGACCACCTGAATGAATTGC
AGTATTTACTCTTCACCATCTTCTTTCTGACCTACATATGCACCTTAGGAGGCAATGTTTTATCATTGT
GGTGACCATAGCTGATTCACCTACACACCCCATGTATTATTCCTAGGAAATCTTGCCCTTATTGAC
ATCTGCTACACTACTACTAATGTCCCCAGATGATGGTGCATCTTCTGTCAGAGAAGAAAATCATTTCCT
ATGGAGGCTGTGTGACCCAGCTTTTGCAATTCATTTCTTTGTTGGCTCAGAGTGTCTCCTCTGGCAGC
AATGGCATATGATCGATATATTGCTATCTGTAAGCCGTTAAGGTACTCATTTATTATGAACAAGGCTCTG
TGCAGCTGGTTAGCAGCTTCATGCTGGACAGGTGGGTTTCTCAACTCAGTGTGACACGGTTTTGACCT
TCCACCTGCCCTTTTGTGGTAACAATCAGATCAATTATTTCTTCTGTGACATACCTCCCTTGCTCATCTT
GTCTTGTTGGTGATACTCCCTCAATGAAGTGGCTTTGCTGTCCATTGGGATCCTCATAGGCTGGACTCCT
TTCCTGTGCGTCATCTTCTACCTTTACATCATCTCCACCATCCTGAGGATCCGTTCTCTGAGGGGA
GGCACAAGCCTTTTCCACCTGTGCCTCCACCTGCTCATTGTTATTCTCTATTATGGCAGTGCTATCTT
CACGTATGTGAGGCCCATCTCATCTTACTCTCTAGAGAAAGATAGATTGATCTCAGTGCTGTATAGTGT
GTCACACCCATGCTGAATCCTGTAATTTATACGCTAAGGAATAAGGACATCAAAGAGGCTGTGAAGGCCA
TAGGGAGAAAGTGGCAGCCACCATTTTCTCTTCTGATATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>MR214842 representing NM_146328 Red=Cloning site Green=Tags(s) MEGKNQTAPSEFIILGFDHLNELQYLLFTIFFLTYICTLGGNVFIIVVTIADSHLHTPMYYFLGNLALID ICYTTTNVPQMMVHLLSEKKIISYGGCVTQLFAFIFFVGSECLLLAAMAYDRYIAICKPLRYSFIMNKAL CSWLAASCWTGGFLNSVLHTVLTFLHPFCGNNQINYYFCDIPPLLILSCGDTSLNELALLSIGILIGWTP FLCVILSYLYIISTILRISSEGRHKAFSTCASHLLIVILYYGSAIFTYVRPISSYSLEKDRLISVLYSV VTPMLNPVIYTLRNKDIKEAVKAIGRKWQPPVFSSDI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_146328
ORF Size:	951 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_146328.2, NP_666440.2</u>
RefSeq Size:	1106 bp
RefSeq ORF:	954 bp
Locus ID:	258325
Cytogenetics:	17 B1
MW:	35.8 kDa