

Product datasheet for MR214858

Cypt8 (NM_001039941) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cypt8 (NM_001039941) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cypt8
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR214858 representing NM_001039941
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTCAGATGGCCAAGAAAGTCCACTGGTCTAGCGCAGCAGCAGGAGCAGCCGCTGCTGCCAAAATTT
 CGAAGCTTGAGAAGACCACAAAAGATTCAAACTATTAAGAAGCGAAACCCTAGTTCTAAGCTTCCCAA
 GAGATCTTCACTCTTTACTTTGTTCTTTCTCGTTCCTGTTGTTGCTGTCGCTGTCGTTGTTGCTGT
 TACTGTCTGTTGCTGCTGTTGTTGCTGTAGTCGTTCTCGTCGTTTTCGTAGCAGAATACTAAGGGTGA
 GATCCCAAGTTCTTTTCAGATTACGGAGAAGGGGGAGCAATCACTTCAAAGAAGAAATTAGGAGACAAT
 GACGCGGAGCCAACTGGAGCTGATCGAACCAGGAAACCGACCATGGCTTTGGAGCCAAAGCGAGATTACAGT
 GCATTCTTCTCATGAGAATGCCAATGTGTCTGAGCCAGAGGAAGTTCCACCATGCCTTGACAGAGACC
 TGGCTCCGAATGGAGAATTGGCCAACTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR214858 representing NM_001039941
 Red=Cloning site Green=Tags(s)

MAQMAKKVHSSAAAGAAAAAISKLEKTTKRFKLIKRNPPSKLPKRSSHLLCSLSRSCCCRCRCCC
 YRCCRCSSRSRRFRSRTLLRVRDPKFFQITEKGEQSLQRRIRRLQTRSQLELIEPEPTMALEPSEITV
 AFFSHENANVSEPEEVPCLDRDLAPNGELANS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9056_b10.zip

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001039941

ORF Size: 519 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039941.2](#), [NP_001035030.1](#)

RefSeq Size: 597 bp


RefSeq ORF: 522 bp

Locus ID: 664723

MW: 19.6 kDa

Gene Summary: DISCONTINUED: This Gene record was discontinued due to lack of supporting evidence.

Product images:



Circular map for MR214858