

Product datasheet for **MR214864**

Fam69c (NM_173770) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Fam69c
Synonyms:	B230399E16Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR214864 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGTGGCCGGGAGGTCAAGAACACGCTAGGCTTGGAGCTGCCAACAACAGCATAGCCCCGCTATGGC
CTGCCAGGCAGGGCCCTGGCTGGCGGCAACAGCTGGCCAGTGCCTGCTCAACAGGAGGAATA
CGTATACTTCAGCCTGCTTCCAGATCTGAGCCGCCACATCCTACCGGTGCTGGGATCCTGCGGTCACTTT
TATGCAGTGGAGTACTTGGCTGCTGGTAGCCCTCACCACAAGGCTCTTCCCACTGGATGATGCTGGCC
AAGCACAGGCTATTAGCCACATTGCCCTCAGCTTCTGGACATGGTGAGCCATTCGACAGTGACTTTTC
CCACCGGCTCCACCTCTGTGATGTCAAACCTGAGAACTTCGCCATCAAGAGGGACTTTACGGTGGTAGCT
ATTGACGTAGACATGGCCTTTTTTGGAGCCTAAATGAGGGAAATCCTTGAGCAGAACTGCACGGGAGATG
AAGATTGCAATTTCTCGATTGTTTCTCAAGTGTGATCTCCGAGTGCACAAGTGTGGAGCCAGCGTGT
CAACAGCAACCTGCAAGTCATCTGTGATAAAATATTCCGACACTGGTTCTCTCAACACACAGGAGCCCG
GCTGTCTCCCTCCAGCTTCGGGTGCAGCTACAGCAGGCCGTGCAGGAGTGTGCTCAGCACGGGGGCTCCT
CGGGGAATTCCTGGACAGCCTCATCCAGTGTGTTCTGGAAGCTTCGATGGCTCTTGCAGGCCACTCTGAA
AGAACTGCAGGAGCGGAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR214864 protein sequence
Red=Cloning site Green=Tags(s)


MVAGEVKNTLGLELPNNSIAPLWPARQGPGRQQLASAWSLQEEYVYFSLPDLSRHILPVLGSCGHF
YAVEYLAAGSPHHKALFPLDDAGQAQAIISHIALSFLDMVSHFSDFSHRLHLCVDPENFAIKRDFTVVA
IDVDMAFFEPKMREILEQNCTGDEDCNFFDCFSKCDLRVHKCGAQRVNSNLQVICDKIFRHWFSSTHRSP
AVSLQLRLQLQQAVQECAQHGGSSGNSWTASSSVFWKLRWLLQATLKELQEAEK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_173770

ORF Size:

792 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min.
	2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
	3. Close the tube and incubate for 10 minutes at room temperature.
	4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
	5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_173770.5</u>
RefSeq Size:	2153 bp
RefSeq ORF:	795 bp
Locus ID:	240479
UniProt ID:	<u>Q8BQT2</u>
Cytogenetics:	18 E4
MW:	29.8 kDa

Product images:



Circular map for MR214864