

Product datasheet for **MR214879L3V**

Bbfl (NM_028377) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Bbfl
Synonyms:	2900006K08Rik; AI427784; Ccdc176
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028377
ORF Size:	1599 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR214879).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_028377.3 , NP_082653.2
RefSeq Size:	2190 bp
RefSeq ORF:	1602 bp
Locus ID:	72873
UniProt ID:	Q3V079
Cytogenetics:	12 D1



Gene Summary:

Basal body protein required in multiciliate cells to align and maintain cilia orientation in response to flow. May act by mediating a maturation step that stabilizes and aligns cilia orientation. Not required to respond to planar cell polarity (PCP) or flow-based orientation cues (By similarity).[UniProtKB/Swiss-Prot Function]