

Product datasheet for **MR214938A1V**

Mouse Olfr1287 (NM_001011773) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr1287 (NM_001011773) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr1287
Synonyms:	MOR248-15
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR214938 representing NM_001011773 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAGAAACAAACCAGACTGTGGTCTCTGAGTTTATTTTCAGGGAATTTGTGCTTCAAAGGAACTAC
AGCTCTTCCTCCTGCTGCCATTTCCATCCTCTATCTAATGGCTGTGGTAGGCAACCTCTTTGTTGTGAT
ATTGATCATCATTGATCATCATCTCCATTCTACCATGTACTTTCTGTTAGCTAATCTATCATTTATTGAC
TTCTGCCTTTCTCAGTTACCACCCCAAACTGATAACTGACCTCCTAAAAGATAACAAAACCATTTCT
TTGGGGCTGCATGAGCCAGATCCTCTGTGTACTTTCTTTGGGGAAGTGAGATGGTGCTTCTTGTAAAC
AATGGCCTATGACCGGTATGTGGCCATCTGCAGGCCACTCCACTACAGCAGCATCATGGGCAGACAGAAG
TGCATCTGGCTTGTGTGATATCATGGACCATTGGCTTTATACATGCCATGAGTCAACTGATACTGATTT
TGGATCTGCCTTTCTGTGGACCTAGAGTGATAGACAGCTTTTCTGTGATATTTCTTTGGTGATGAAATT
AGCCTGCATGAATACTGATACTCTGGAATTCTAATAAATGCTGACAGTGGTGTGTTGGCAACAATTGT
TTCATCCTCTTACTGATATCTTATACTAACATTCTATTAAGTGTCACTTCACTCTAAAGATGGTTCAT
CAAAGGCATTATCTACATGCACATCCCATATCATAGTGGTTTTGCTGTTCTTTGGGCCAGTCATTTTCAT
CTATCTGTGTCCAGTTAGCATCACCTGGGTTGACAACCTTTCTTGCCGATTTTACTCAGTCATCACACCT
CTCCTTAATCCAGCCATCTACACACTGAGAAACAAAGATATTAATAAATGCCATAAAGAAGCTGATAAACC
ACTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR214938 representing NM_001011773 Red=Cloning site Green=Tags(s) MEETNQTVVSEFIFQGICASKELQLFLLLPFSILYLMVVGNLFVVILIIDHHLHSTMYFLLANLSFID FCLSSVTTPKLITDLLKDNKTISFGGCMSEQILCVYFFGGSEMVLLVTMAYDRYVAICRPLHYSSIMGRQK CIWLVISWTIGFIHAMSQILILDLPCGPRVIDSFCDISLVMKLACMNTDTLEILINADSGVLATTC FILLISYTNILLTVQLHSKDGSSKALSTCTSHIIVLLFFGPVIFIYLCPVSITWVDNFLAVFYSVITP LLNPAIYTLRNKDIKNAIKKLINHL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001011773
ORF Size:	915 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001011773.1, NP_001011773.1</u>
RefSeq Size:	918 bp
RefSeq ORF:	918 bp
Locus ID:	257935
Cytogenetics:	2 E3
MW:	34.3 kDa