

Product datasheet for **MR214945**

Olfr1276 (NM_146395) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1276 (NM_146395) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr1276
Synonyms:	MOR245-10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR214945 representing NM_146395 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAATGAAATAAATTACACCGAGGTGTCTGAATTTGTGTTTCTGGGACTTTCAACATCTAAACACATAC
 AGCATTCTTTCTTGCCTTCTCTGTGGTATTTATGTGACCATTGTGCTGGGCAATACTCTTGTGTATT
 TACTCTAGCCTTTGACCCACATTTACATTCCCCTATGTACTTTCTTTGGGAAACCTCTCATTTATTGAC
 TTATGTCTTTCCACCTTAACCGTACCCAAGATGATTTCTGACCTGTCCTCTGGGCATAATACCATCTCAT
 TCCAGGGTTGTGCTTCCAGATATTTGTTCTTCATGTTCTTGGGGCATCTGAGATGGTATTGCTGGTAGC
 CATGGCCTGGGACAGATATGTGGCCATATGCAAAACCCTCCACTATTTAACCATCATGAGCCACGGATG
 TGCTATTGCTTCTATCTGGGGCTTGATTATTGGCTTTTACATTCACTGGCCCACTAGGTTTGTG
 TCCATTGAGATTTTGTGGACCTAATAAGATAGATAGTTTTTATTGTGATCTTCTAGGTTTATCAAAC
 GGCTGCATAGACAACCTACAGAAATGGAGTTCATGGTTGCTGCCAACAGTGGCATCATTTCTATTGGCACC
 TTCTTCTTACTGATAATCTCCTACATTGTTATCCTGTTAATGTAAGAAAACATTCATCAGGAGATTTGT
 CCAAGGCCCTCTCAACACTTTCTGCTCACATCTCCGTGGTAGTTTTGTTCTTTGGACCATGCATCTTCGT
 GTACATGTGGCCGTTTCTACTGTGCCAGTGGATAAGTTCCTTGCCATTCTGGACTTCATGATTACCCC
 ATCCTGAACCTGCCATTTACACGCTGAGGAACAAAGACATGAAGGTGGCAATGAGGAACTGAGTTATC
 AGTTCCTGAATTTAGGAAAATGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR214945 representing NM_146395
 Red=Cloning site Green=Tags(s)

MNEINYTEVSEFVFLGLSTSKHIQHFFLAFSVVFYVTIVLGNLTVFTLAFDPLHSPMYFLLGNLSFID
 LCLSTLTVPKMISDLSSGHNTISFQGCVFQIFVLHVLGASEMVLLVAMAWDRYVAICKPLHYLTIMSPRM
 CLLLLSGAWIIGFLHSAQLGFVVHLRFCGPNKIDSFYCDLPRFIKLACIDNYRMEFMVAANSGLISIGT
 FFLLIISYIVILFNVRKHSSGDLKALSTLSAHISVVVLFPGPCIFVYMWPFTVPVDKFLAILDFMITP
 ILNPAIYTLRNKDMKVMRKLSYQFLNFRKMS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9067_b02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_146395

ORF Size: 936 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_146395.1](#), [NP_666507.1](#)

RefSeq Size: 939 bp

RefSeq ORF: 939 bp


Locus ID: 258390

Cytogenetics: 2 E3

MW: 35.4 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR214945