

Product datasheet for MR214961

Spin2-ps1 (NM_029106) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spin2-ps1 (NM_029106) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Spin2-ps1
Synonyms:	Gm14354; OTTMUSG00000016414
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR214961 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGAAAGTTCAAAGCTGACCAAGAAAAGAGCTGGTCGGAAGAGGCAGAGGAGCAGATCTCCAGCCCTGC
 CCAAAGGAACATCTTGGGCTGCAGATTTCTCATAAGTGGAAGAAGGAGATGAGTCTATCACACAGTG
 GAATGGAAGTGTGCTGGATCAGGTGCCAGTAAACCCTTCCTGTACCTGGTTAAATATGATGAGATTGAC
 GCTGTGCATGCTCTGGAACCTTACAAAGATAAGAGGGTATTGTCCCTAAAAGTTATCGCCAAGAGGGTGG
 TATCATCCGGAGTCACTGATTCCAGCTTTGTTGATGCCATAGTGGCAAAGAAGTTAATCATCTATTTGA
 GGGTGAGCATGGCTCTAAGGAAGAGTGGAGGGGGATGGTCCTGAGCCAAGCGCCTATCTTGATAACAGC
 TTTTATATTACTTATGAGAGGGATCCCATTATACACGTACCAGCTTCTGGATGACTTTAAAGAGGGTG
 ACCTACAAATCACGGAGGGGATCAGTGACCCCCCTTCTTTAGACATCGACCTGGAGCTTGTGGATGGCCT
 GATAGGAAAACACGTGGAACACAAATGATGACGGCTCAAGAGGGATGGCCTGATCATTTACCAGATT
 GAAACCAACCTAGAGTGTACCTCATCAAGTATGAAGATGATGTTTATATACATGTAACCTCACTTGGAGA
 AAGAGTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)


Protein Sequence: >MR214961 protein sequence
 Red=Cloning site Green=Tags(s)

MESSKLTKKRAGRKRQSRSPALPKRNILGCRISHKWKEGDESITQWNGTVLDQVPVNPSTLYLVKYDEID
 AVHALEYKDKRVL SLKVI AKRVVSSGVTSSSFVDAIVGKEVNHLFEGEHGSKEEWRGMVLSQAPILDNS
 FYITYERDPILYTYQLDDFKEGDLQITEGIDPPSLDIDLELVDGLIGKHVENTNDDGSKRDGLIYQI
 ETKPRVYLIKVEDDVHIVHLEKEF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_029106

ORF Size: 708 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_029106.1](#), [NP_083382.1](#)

RefSeq Size: 778 bp


RefSeq ORF: 711 bp

Locus ID: 74851

Cytogenetics: X A1.1

MW: 27 kDa

Product images:



Circular map for MR214961