

Product datasheet for **MR214974L3V**

Syce2 (NM_027954) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Syce2 (NM_027954) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Syce2
Synonyms:	1700013H19Rik; AA407907; Cesc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027954
ORF Size:	513 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214974).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027954.3 , NP_082230.2
RefSeq Size:	911 bp
RefSeq ORF:	516 bp
Locus ID:	71846
UniProt ID:	Q505B8
Cytogenetics:	8 C3


[View online »](#)

Gene Summary:

Major component of the transverse central element of synaptonemal complexes (SCS), formed between homologous chromosomes during meiotic prophase. Requires SYCP1 in order to be incorporated into the central element. May have a role in the synaptonemal complex assembly, stabilization and recombination.[UniProtKB/Swiss-Prot Function]