

Product datasheet for MR214988

Tmem40 (NM_001168257) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem40 (NM_001168257) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tmem40
Synonyms: 9030407H20Rik; AV091463; W64236
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR214988 representing NM_001168257
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGCTTCAGGATCTTCTCCAGTCTCAAGACAGTGGTGGAGTCCACAGGAAACGGAAGATCACT
 ACCAAGAGACAGAGCTCCACAAGCATCATGGGAAGGCCCGGAGCGATATAAGAGAGACAAGTCTTCCTC
 TTCATCCTCATCCTCGTCATCTTCTCCTCTCTCTCCTCATCTTCTCCTCCTCAGGTGCCGAC
 CATGGGAACTGGAAGTCTAAAGGACGAGCTCAGCTCTGCGGAGGTGCTGCTGGGAAATGGTGCCTA
 CTGGGAGTCAGGACTCCGAAGGAGAGGTTCCGGCTCAGCGAGGAGAGAAGTGGAGGCCTCTCAGTTAAG
 AAGACTGAACATAAAGAAGGATGATGAGTTTTTCCATTTTGCCTTCTCTGCTTTGCCATCGGAGCCTTG
 CTGGTGTGTTACCACTACTATGCAGACTGGTTCATGTCCCTCGGGGTCCGGCTGCTCACCTTTGCTTCTC
 TCGAGACCATTGGCATCTATTTCCGGCTCGTGTACCGGATCCACAGCGTGCTTCAGGGCTTCATCCCCCT
 CCTTCAGAAGTCCGGCTGCCAGGGTTCAGGAGAACTAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR214988 representing NM_001168257
 Red=Cloning site Green=Tags(s)

MEASGSSSQSDSGGVHRETEDHYQETELHKHHGKARERYKRDKSSSSSSSSSSSSSSSSSSSSSSGAD
 HGELEVLKDELQLCGGAAGEMVPTGESGLRRRSGSAEVEASQLRRLNIKKDDEFFHFVLLCFAIGAL
 LVCYHYADWFMSLVGLLTFASLETIGIYFGLVYRIHSVLQGFIPLLQKFRPLPGFRRTN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001168257

ORF Size: 600 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001168257.1](#), [NP_001161729.1](#)

RefSeq Size: 1339 bp


RefSeq ORF: 603 bp

Locus ID: 94346

Cytogenetics: 6 E3

MW: 22.3 kDa

Product images:



Circular map for MR214988