

## Product datasheet for **MR214994L3V**

### 1700047117Rik2 (NM\_001100116) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	1700047117Rik2 (NM_001100116) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	1700047117Rik2
Synonyms:	MGC61221
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001100116
ORF Size:	621 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR214994).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001100116.1</a> , <a href="#">NP_001093586.1</a>
RefSeq Size:	3863 bp
RefSeq ORF:	624 bp
Locus ID:	100101807
UniProt ID:	<a href="#">Q8BR63</a>
Cytogenetics:	12 C1



[View online »](#)