

Product datasheet for **MR215014L4V**

1700028K03Rik (NM_182745) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	1700028K03Rik (NM_182745) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	1700028K03Rik
Synonyms:	MGC117884
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182745
ORF Size:	558 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215014).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182745.2 , NP_877422.1
RefSeq Size:	795 bp
RefSeq ORF:	561 bp
Locus ID:	76421
UniProt ID:	Q3KQP7
Cytogenetics:	5 F



[View online »](#)

Gene Summary:

Plays a key role in reinforcing the integrity of the central element of the synaptonemal complex (SC) thereby stabilizing SC, ensuring progression of meiotic prophase I in male and female germ cells (PubMed:30949703). Promotes homologous recombination and crossing-over in meiotic prophase I via its association with SHOC1 (PubMed:30746471). Required for the localization of TEX11 and MSH4 to recombination intermediates (PubMed:30746471). [UniProtKB/Swiss-Prot Function]