

Product datasheet for **MR215041**

Yae1d1 (NM_025904) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Yae1d1
Synonyms:	1600012F09Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR215041 representing NM_025904
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCTTGGTTTCGAGCGGCTCCGGAGAACGCGAGCCCTGGAGACCAGGACGTGTTGACGAGGAAGCGG
ACGAGTCTCTGTTAGCGCAGCGGGAATGGCAGGGTCACATGCGGAAACGAGTCCAAGAGGGTTACAGAGA
TGGGCTAGATGCCGGCAAAGCACTTACTCTTCAACGGGGCTTCAATCAAGGGTACAAAGAGGGTGCAGAC
GTCATTATAAACTATGGACTCCTTAGAGGAACACTGAGTGCTTTGCTTTCTGGTGTACCTTCATGGTA
ATAGTTCGACTCTGATCAGTAAGATAAGTAACCTCCTGGACACAGTTGACCAGTGTGAAGAGTGTGTGCT
CAGACGTCTGAAATCGGTTACATCACAGCCACGAGTTACGGATTTATTGGACGCCATTGAGGATATGGAC
CTTTGTCACGTAGTTCGGTCAGACCAAGCGAAAGGTGACAGTAACGCTGAGATTAACACCGATTCTGACA
GTCCAAATGGGATGGGTTGTTACATCCAGAACATTGTAGACCAAAGAGCATGCACATTTGAAAAAACC
AAGCCTCTCTTGATTCTGGAGCAGACATCCAGTTAGTGAACAGCTGGGAGTACCGGCAGACATATTA
CAGCACCTGAGACAACTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR215041 representing NM_025904

Red=Cloning site Green=Tags(s)

MSWFRAAPENASPGDQDVFDDEADESLLAQREWQHMRKRVQEGYRDGLDAGKALTQRFNQGYKEGAD
VIINYGLLRGTL SALLSWCHLHGNSSTLISKISNLLDTVDQCEECVLRRLKSVTSQPRVTDLLDAIEDMD
LCHVVPDQAKGDSNAEINTDSDSPNGMGCSHPEHCRPKEHAHLEKPSLSWILEQTSSLVKQLGVPADIL
QHLRQL

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mm9059_f04.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_025904

ORF Size:

648 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_025904.3](#), [NP_080180.1](#)

RefSeq Size: 3436 bp


RefSeq ORF: 651 bp

Locus ID: 67008

Cytogenetics: 13 A2

MW: 24 kDa

Product images:



Circular map for MR215041