

Product datasheet for **MR215041L3V**

YaeId1 (NM_025904) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	YaeId1
Synonyms:	1600012F09Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025904
ORF Size:	648 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR215041).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_025904.3](#), [NP_080180.1](#)

RefSeq Size: 3436 bp

RefSeq ORF: 651 bp

Locus ID: 67008

Cytogenetics: 13 A2

