

## Product datasheet for MR215079

### Olfr1065 (NM\_146408) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1065 (NM_146408) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr1065
Synonyms:	MOR190-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR215079 representing NM_146408 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGAACAAAACCGCACAGTAGTAACTGAATTCATCCTGATGGGCATCACTGACCACCCTGAGCTGC  
AATCCTCATTGTTTGGTTTGTTCCTTCATCTATCTGATCACTGGTGGGAAATTTGGGCATGATTGT  
CCTTACAATGGTGGACTCCAGGCTGCAACACCCATGATTTCTTCTTAGACACCTCGCTACTACTGAC  
CTTGTTATTCAACAACGTGGGACCCAAAATGTTAAACAATTTATTGTTGATCAAAAATAGAATATCCT  
TTAACCTTTGTGCTACACAGTTAGCCTTCTTCTCTGTTTCATTGCTTGTGAACCTTTTATTCTGTCTGC  
AATGTCCTATGACCGCTATGTGGCTATCTGTAAGCCTCTGTTCTATATGGTCATCATGTCAAAAAGACTA  
TGCTGGGTGCTTGTAGTAATCCCTTATGTCTACTGTACAATTGTGGCTCTTTAATCACAATTAAGATT  
TCACTTTGTCTTTCTGTGGCTCCAATGTCATCAGTCATTTCTACTGTGACAGTCTTCCTTTGTTATCTTT  
GGTCTGCTCAAAATACACAAGAAATTGAAGTCATAATGTTGTTCTTATCAGCTTTTAAATTTAATTCCTCT  
CTTTTGCTTGTCTTGTGTCCTACCTGCTCATCCTTATAGCTATACTTAGGATGAACTCTGCTGAGGGCA  
GGCGCAAGGCTTCTCCACCTGTGGATCCCACCTGACTGTGGTCACTGTCTTCTATGGAACGTTGATATT  
TATGTACATGCAGCCAAATCCAGTCATTCCTTTGACACTGATAAATTAGCTTCTGTTTTTACACCCTG  
ATTATCCCATGTTGAATCCCTTGATCTACAGCTTGAGGAATAAAGATGTAAAAGATGCTCTTCATAGGA  
CAGGGAAAAAACTACACAGTGTGGCTTT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR215079 representing NM\_146408  
Red=Cloning site Green=Tags(s)

MENQNRVTVEFILMGITDHPQLQSSLFGLFLLIYLISLVGNLGMIVLTMVDSRLQTPMYFFLRHLATTD  
 LGYSTTVGPKMLTNFIVDQNRISFNLCATQLAFFLLFIACELFILLSAMSYDRYVAICKPLFYMVIMSKRL  
 CWVLVVIPIVYVYCTIVALLITIKIFTLFCGNSVISHFYCDLPLLSLVCSNTQEIEVIMLFLSAFNLISS  
 LLLVLSVYLLILIAILRMNSAEGRRKAFSTCGSHLTVVTVFYGTILFMYMQPKSSHSFDTDKLASVFTL  
 IIPMLNPLIYSLRNKDVKDALHRTGKKLHSVCF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9067\\_b08.zip](https://cdn.origene.com/chromatograms/mm9067_b08.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_146408

**ORF Size:** 939 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_146408.2](#), [NP\\_666520.2](#)

**RefSeq Size:** 992 bp

**RefSeq ORF:** 942 bp


**Locus ID:** 258403

**Cytogenetics:** 2 D

**MW:** 35.5 kDa

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

**Product images:**



Circular map for MR215079