

Product datasheet for MR215125

Wdr33 (NM_001170966) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Wdr33 (NM_001170966) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Wdr33
Synonyms: 1110001N06Rik; 2310011G05Rik; 2810021O11Rik; 8430413N20Rik; WDC146
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR215125 representing NM_001170966
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTACAGAAATTGGCTCTCCTCCACGATTTTTCCATATGCCGAGTTTCAACATCAGGCTCCCCGAC
AATTATTTTATAAGCGACCTGATTTTGCACAGCAGCAAGCAATGCAGCAGCTCACCTTTGATGGAAAACG
AATGAGAAAAGCTGTAAACCGAAAACCATAGACTACAATCCCTCTGTAATTAAGTATTTGGAGAATAGG
ATATGGCAAAGAGATCAGAGAGATATGCGGGCAATTCAGCCGGATGCAGGTTATTATAATGATCTGGTCC
CACCAATAGGAATGTTGAATAATCCTATGAATGCAGTAACAACAAAATTTGTTAGGACATCAACAAACAA
AGTAAATGCCAGTATTTGTTGTGAGGTGGACCCCGAAGGAAGACGATTGGTCACTGGAGCTTCTAGT
GGGGAGTTCACCTTGTTGAATGGGCTGACTTTCAACTTTGAAACCATATTGCAGGCTCATGATAGCCCTG
TGAGGGCCATGACTTGGTCACATAATGACATGTGGATGTTGACAGCAGATCATGGAGGATATGTGAAATA
TTGGCAATCCAACATGAACAACGTCAAGATGTTCCAGGCACATAAGGAGGCGATTAGAGAGGCCAGGTTT
ATACACAATATACCATTTTCTGTAGTCCCTATTGTCATGGTTAAATTAATCTCTAAGTGTATTCTGGGTG
CCAAGATGCATGGGCTCTCTCAGATTCTGGGAACTTTCTGCACCTTATGAACACAATTTCTCTTTGT
TTTCACACATTCACATTTTCTGGCACCTTTCTGAAGTAGTTTTGTCCTGGTATCAGCCTTTGCAATAT
GTTAGAGATGATTGTCTGCCACATTTTGCACCTGGTTTTCTCTTTTATTATGGCTAATAATGTGTATAC
GTTATTCTTTTTATTATCTACTGTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001170966.1 , NP_001164437.1
RefSeq Size:	2548 bp
RefSeq ORF:	939 bp
Locus ID:	74320
Cytogenetics:	18 B1
MW:	37 kDa
Gene Summary:	Essential for both cleavage and polyadenylation of pre-mRNA 3' ends.[UniProtKB/Swiss-Prot Function]