

Product datasheet for MR215126

Wdr33 (NM_001170967) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Wdr33 (NM_001170967) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Wdr33
Synonyms: 1110001N06Rik; 2310011G05Rik; 2810021O11Rik; 8430413N20Rik; WDC146
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR215126 representing NM_001170967
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTACAGAAATTGGCTCTCCTCCACGATTTTCCATATGCCGAGTTTCAACATCAGGCTCCCCGAC
AATTATTTTATAAGCGACCTGATTTTGCACAGCAGCAAGCAATGCAGCAGCTCACCTTTGATGGAAAACG
AATGAGAAAAGCTGTAAACCGAAAACCATAGACTACAATCCCTCTGTAATTAAGTATTTGGAGAATAGG
ATATGGCAAAGAGATCAGAGAGATATGCCGGCAATTCAGCCGGATGCAGGTTATTATAATGATCTGGTCC
CACCAATAGGAATGTTGAATAATCCTATGAATGCAGTAACAACAAAATTTGTTAGGACATCAACAAACAA
AGTAAATGCCCAGTATTTGTTGTGAGGTGGACCCGGAAGGAAGACGATTGGTCACTGGAGCTTCTAGT
GGGGAGTTCACCTTGTTGAATGGGCTGACTTTCACTTTGAAACCATATTGCAGGCTCATGATAGCCCTG
TGAGGGCCATGACTTGGTCACATAATGACATGTGGATGTTGACAGCAGATCATGGAGGATATGTGAAATA
TTGGCAATCCAACATGAACAACGTCAGATGTTCCAGGCACATAAGGAGGCGATTAGAGAGGCCAGTTTC
TCACCCACGGATAATAAATTTGCTACATGCTCTGATGACGGCACTGTTAGAATCTGGGATTTTCTTCGTT
GCCATGAGGAAAGAATTCTCCGAGGGTTTCTTTGTAGCCCTGGCTGTCTGAAATTTGCTCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR215126 representing NM_001170967
 Red=Cloning site Green=Tags(s)

MATEIGSPRRFFHMPRFQHQAPRQLFYKRPDFAQQQAMQQLTFDGKMRKAVNRKIDYNPSVIKYLENR
 IWQRDQRDMRAIQPDAGYYNDLVPIGMLNPNMNAVTTKFKVVRTSTNKVKCPVFFVVRWTPGRRLLVTGASS
 GEFTLWNGLTFNFETILQAHDSPVRAMTWSHNDMWMLTADHGGYVKYQSNMNNVKMFQAHKEAIREASF
 SPTDNKFATCSDDGTVRIWDFLRCHERILRGFLFVALAVLEIAQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001170967
 ORF Size: 765 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001170967.1 , NP_001164438.1
RefSeq Size:	1148 bp
RefSeq ORF:	768 bp
Locus ID:	74320
Cytogenetics:	18 B1
MW:	30.2 kDa
Gene Summary:	Essential for both cleavage and polyadenylation of pre-mRNA 3' ends.[UniProtKB/Swiss-Prot Function]