

Product datasheet for **MR215134L3V**

Abi3bp (NM_001014399) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Abi3bp (NM_001014399) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abi3bp
Synonyms:	5033411B22Rik; AI506287; BG172926; D930038M13Rik; eratin; TARSH
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001014399
ORF Size:	2727 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215134).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001014399.2 , NP_001014421.2
RefSeq Size:	3884 bp
RefSeq ORF:	2730 bp
Locus ID:	320712
Cytogenetics:	16 C1.1



[View online »](#)