

Product datasheet for **MR215153**

Slc24a2 (NM_172426) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc24a2 (NM_172426) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc24a2
Synonyms:	2810021B17Rik; 6330417K15Rik; AI847460; Nckx2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR215153 representing NM_172426
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATCTCCACCAAAGCCCCACAGCCAGGCTCCTGCAGAAATGGTGTTCCCATGAGTCACCATTGGCT
 GCAGGAGACATTATAATAGCAGGAAAAAACTGAAGCTAATTCGAGTCATTGGCCTTGTCATGGGCCTGGT
 AGCTGTGAGCACTGTCCATTTTCAATCAGTGCCTTTACTGAGACAGACTCTCAGAGCAACAGAGGAGAA
 GCCAGTGATATGAGTGGCCCCAGAGTAGCACAGGGTCACCGTCAAAGAACTCTCTTAGACTTAAATGACA
 AGATTCGAGATTACCCCCACAGCCACCTGCTTCTCAGGAAGACCAAGCTGAGAACAGTACTGAGCACAC
 CCAGGGGAGACTACCCAAAGGACATCTTTCCCTGGAGGAGCGAAGGAAAGGTGCCATCATTCTGCATGTC
 ATTGGGATGATCTATATGTTTCATAGCATTAGCCATCGTCTGTGATGAGTTCTTTGTTCTCTCTAACTG
 TCATCACTGAAAACTGGGAATCTCTGATGATGTGGCTGGAGCCACCTTCATGGCTGCAGGCGGGTCAGC
 CCCAGAACTTTTACGTCTCTTATAGGGGTCTTCATTGCTCACAGCAATGTCGGCATTGGCACCATCGTG
 GGTTCCGCCGTGTTCAATATTCTTTGTTATTGGCATGTGTGCTCTGTTTTCTAGAGAAATCCTAAACC
 TGACCTGGTGGCCGCTCTTCGGGATGTGTCCTTCTACATCGTTGATTGCTCATGTTGATCACATTTTT
 CCTGGATAACGTTATCATGTGGTGGGAAAGCTTACTCCTCTTAACAGCTTATTTGCGCTATGTGGTTTTT
 ATGAAGTTCAACGTCCAAGTAGAGAGATGGGTGAAGCAATGATATCTCGCAATAATGTTATCAAGGTGA
 CAGTGCCAGAAGCCCAAGCAAAGTACCTACAGCTGGGGACAAGGATGGGCCAACTTTACCGTCGAAGCC
 ACGCCTCCAACGAGGCGGAAGCTCTGCCTCCCTCCACAACAGTCTCATGAGGAATAGCATCTTCCAGCTC
 ATGATACACACCTTGACCCACTCGCTGAAGAAGTTGGATCATATGAAAACTAAATATTATGACACAA
 TGACTGAAGAAGGAGGTTTCAGAGAAAAAGCTTCAATTCTCCACAAGATCGCCAAGAAGAAATGCCAAGT
 GGATGAGAATGAGAGGCAGAATGGGGCTGCCAATCATGTGAAAAAATTGAACTCCCAACAGCACAAAGC
 ACAGAGGTGGAATGACGCCATCGAGCGAAGCTTCGGAACCCGTCGAGAATGGAACCTCTCCACAACA
 TTGAAGCTGCAGATGCCCGAAGGCCACAGAGACTGCCGAGGAAGAAGATGACCAGCCCCCTCAGCTTGTC
 CTGGCCGACCAACACCCGCAAGCAAGCCACCTTCTAATCGTTTTTCCATTGTGTTCCCTCTCTGGATC
 ACCTTACCTGACGTCCGAAAACCTGCTTCCAGGAAGTTCTTCCCTATCACATTCTTCGGCTCCATCACCT
 GGATTGCAGTTTTCTCATACCTGATGGTGTGGTGGGCACACCAGGTTGGAGAGACAATCGGCATTAGTGA
 AGAAATCATGGGTCTGACCATCTTGGCTGCTGGGACCTCTATCCCTGATCTTATCACCAGTGTATAGTG
 GCCCGAAAGGGCCTGGGGACATGGCTGTGTCCAGCTCTGTTGGAAGCAACATCTTTGACATCACTGTGG
 GGCTCCCACTGCCCTGGCTCCTCTACACCATCATCCACAGATTACGCCAGTGACTGTGACGAGCAATGG
 CCTCTTCTGTGCCATCGTCCTTCTTTCATCATGTTGCTTTTGTCACTTCTCCATCGCCCTCTGCAAA
 TGGCGGATGAACAAAATCCTCGGCTTATCATGTTTGGCCTCTACTTCGTGTTCTGGTGGTCAGCGTCC
 TCCTGGAAGACAAAGTTCTCGTGTGCCCTGTCTCCATC

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215153 representing NM_172426
Red=Cloning site Green=Tags(s)


MDLHQSPtarLLQKWCShESpFGcRRhYNSrKKLkLIRvIGLVMGLVAVSTVPFSISaFTETDSQSNRGE
ASDMSGPRVAQgHRQRTLLDLNDKIRDYTPPPASQEDQAENSTEHTQGDYPKDIFSLERrKGaIILHV
IGMIYMFIALaIVCDEFFVPSLTVITEKLGISDDVAGATfMAAGGSaPELFTSLIGVfIAHSNVGIGTIV
GSAVFNILfVIGMCALfSREILNLtWWPLFRDVSfYIVDLLMLITffLDNVIMWWESLLLLtAYfAYVfV
MKFNvQVERWVKMISrNNVIKVTvPEaQAKSPTAGDKDGPTLSPKPrLQRGGSSaSLHNSLMrNSIfQL
MIHTLDPLaEELGSYGKLKYyDTMTEEGrFREKASILhKIAKKKQVDENERQNGaANHVEKIELPNSTs
TEVEMTPSSEaEPVQNGNLShNIEaADAPKaTETAEEEDDQPLSLSWPTNTRKQATfLIVfPIVfPLW
TLpDVRKPaSRKFFPITffGSITWIAVfSYLMVWwAHQVGETIGISEEIMGLTILAAGTSIPDLITSVIV
ARKGLGDMAVSSSVGSNIfDITVGLPLPwLLYtIIHRfSPVTvSSNGLfCAIVLLfIMLLfVILSIALCK
WRMNkILGFIMfGLYfVfLVVSVLLLEDKVLVCPVSI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9004_h08.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_172426

ORF Size: 1998 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_172426.2](#), [NP_766014.1](#)

RefSeq Size: 10606 bp


RefSeq ORF: 2001 bp

Locus ID: 76376

Cytogenetics: 4 C4

MW: 74.7 kDa

Product images:



Circular map for MR215153