

## Product datasheet for MR215175L4

### Kng2 (NM\_001102410) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Kng2 (NM_001102410) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Kng2
Synonyms:	MGC117836; MGC123404
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215175).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 385643

Cytogenetics: 16 B1