

## Product datasheet for **MR215182L3V**

### Ank3 (NM\_009670) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ank3 (NM_009670) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ank3
Synonyms:	2900054D09Rik; AI314020; An; Ank; Ank-3; AnkG; Anky; Ankyrin-3; Ankyrin-G
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009670
ORF Size:	2694 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215182).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009670.4</a> , <a href="#">NP_033800.2</a>
RefSeq Size:	6966 bp
RefSeq ORF:	2697 bp
Locus ID:	11735
UniProt ID:	<a href="#">G5E8K5</a>
Cytogenetics:	10 36.1 cM



[View online »](#)

**Gene Summary:**

This gene encodes a member of the ankyrin protein family. Ankyrins link integral membrane proteins to the spectrin-based cytoskeleton. Ankyrin family members share a protein structure which includes three independently folded domains: the N-terminal ankyrin repeat domain, the central spectrin-binding domain, and the C-terminal rod domain. This ankyrin functions as the major ankyrin in the kidney and may play a role in the polarized distribution of many integral membrane proteins to specific subcellular sites. Alternative splicing of this gene results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]