

## Product datasheet for **MR215216L4V**

### 1110017D15Rik (NM\_001048005) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	1110017D15Rik (NM_001048005) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	1110017D15Rik
Synonyms:	Cbe1; Smrp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001048005
ORF Size:	780 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215216).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001048005.1</a> , <a href="#">NP_001041470.1</a>
RefSeq Size:	1127 bp
RefSeq ORF:	783 bp
Locus ID:	73721
UniProt ID:	<a href="#">Q2MH31</a>
Cytogenetics:	4 A5
Gene Summary:	May be involved in differentiation or function of ciliated cells (By similarity). May play a role in spermatogenesis.[UniProtKB/Swiss-Prot Function]



[View online »](#)