

Product datasheet for **MR215240L4V**

Fancf (NM_001115087) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fancf (NM_001115087) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fancf
Synonyms:	A730016A17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001115087
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215240).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001115087.1 , NP_001108559.1
RefSeq Size:	1691 bp
RefSeq ORF:	1032 bp
Locus ID:	100040608
UniProt ID:	E9Q5Z5
Cytogenetics:	7 B4
Gene Summary:	DNA repair protein that may operate in a postreplication repair or a cell cycle checkpoint function. May be implicated in interstrand DNA cross-link repair and in the maintenance of normal chromosome stability.[UniProtKB/Swiss-Prot Function]



[View online »](#)