

Product datasheet for MR215331

Sprr2a1 (NM_011468) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sprr2a1 (NM_011468) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Sprr2a1
Synonyms: Sprr2a; Sprr2a3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR215331 representing NM_011468
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTTACTACCAGCAGCAGTGAATCAGCCGTGCCGGCCTCCTCTGTGTGCCACCCCGAAGTGCC
CTGAGCCTTGTCTCCCAAGTGTGGCCTGGGCCTTGTGTCCTGTATGTGCTTTGAGCCTTGTCTTCC
TTCAGTGTGGCCTGGCCGTGTCGTCCTGTAGTGTGCTATGAGCAATGCCACCTCAGCCATGGCAGTCA
ACATGCCCTCCTGTGCAATTTCCACCATGCCAGCAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215331 representing NM_011468
Red=Cloning site Green=Tags(s)
MSYYQQCNQPCRPVPCPPKCEPCPPQVWPGPCRPVMCFEPCLPSVWPGPCRPVVCYEQPPQPWQS
TCPVQFPQCQK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9076_c05.zip

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_011468

ORF Size: 249 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011468.4](#), [NP_035598.2](#)

RefSeq Size: 3547 bp

RefSeq ORF: 252 bp

Locus ID: 20755


UniProt ID: [Q4KL71](#)

Cytogenetics: 3 40.14 cM

MW: 9.4 kDa

Gene Summary: Cross-linked envelope protein of keratinocytes. It is a keratinocyte protein that first appears in the cell cytosol, but ultimately becomes cross-linked to membrane proteins by transglutaminase. All that results in the formation of an insoluble envelope beneath the plasma membrane (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR215331